
|                 |                                                                     |
|-----------------|---------------------------------------------------------------------|
| <b>English</b>  | User Manual <b>Analytical and Precision Balances</b> MR             |
| <b>Deutsch</b>  | Benutzerhandbuch <b>Analysen- und Präzisionswaagen</b> MR           |
| <b>Español</b>  | Manual de usuario <b>Balanzas analíticas y de precisión</b> MR      |
| <b>Français</b> | Guide de l'utilisateur <b>Balances d'analyse et de précision</b> MR |
| <b>简体中文</b>     | 用户手册 <b>分析与精密天平</b> MR                                              |
| <b>日本語</b>      | ユーザマニュアル <b>分析天びんと上皿天びん</b> MR                                      |
| <b>Čeština</b>  | Návod k použití <b>Analytické a přesné váhy</b> MR                  |
| <b>Dansk</b>    | Brugervejledning <b>Analyse- og præcisionsvægte</b> MR              |
| <b>Hrvatski</b> | Korisnički priručnik <b>Analitičke i precizne vage</b> MR           |
| <b>Magyar</b>   | Felhasználói útmutató <b>Analitikai és precíziós mérlegek</b> MR    |

**METTLER TOLEDO**

► [www.mt.com/MR-RM](http://www.mt.com/MR-RM)

en For full information, always refer to the Reference Manual (RM).

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# Table of Contents

|          |                                                         |           |
|----------|---------------------------------------------------------|-----------|
| <b>1</b> | <b>Introduction</b>                                     | <b>3</b>  |
| 1.1      | Document purpose .....                                  | 3         |
| 1.2      | Further documents and information .....                 | 3         |
| 1.3      | Acronyms and abbreviations .....                        | 3         |
| <b>2</b> | <b>Safety Information</b>                               | <b>4</b>  |
| 2.1      | Definition of signal words and warning symbols .....    | 4         |
| 2.2      | Product-specific safety notes .....                     | 4         |
| <b>3</b> | <b>Design and Function</b>                              | <b>5</b>  |
| 3.1      | Overview analytical balances .....                      | 5         |
| 3.2      | Overview precision balances, with draft shield .....    | 6         |
| 3.3      | Overview precision balances, without draft shield ..... | 7         |
| 3.4      | Overview terminal .....                                 | 7         |
| 3.5      | Overview interface connections .....                    | 8         |
| 3.6      | Components description .....                            | 8         |
| 3.6.1    | Draft shield .....                                      | 8         |
| 3.6.2    | Weighing pan .....                                      | 9         |
| 3.6.3    | Drip tray .....                                         | 9         |
| 3.6.4    | Door handle .....                                       | 9         |
| 3.6.5    | Leveling feet .....                                     | 9         |
| 3.6.6    | Terminal .....                                          | 10        |
| 3.6.7    | QuickLock for draft shield .....                        | 10        |
| 3.6.8    | QuickLock for doors and front panel .....               | 10        |
| 3.6.9    | Release button for back panel .....                     | 10        |
| 3.7      | User interface .....                                    | 11        |
| 3.7.1    | Main sections at a glance .....                         | 11        |
| 3.7.2    | Main weighing screen .....                              | 12        |
| <b>4</b> | <b>Installation and Putting into Operation</b>          | <b>13</b> |
| 4.1      | Selecting the location .....                            | 13        |
| 4.2      | Unpacking the balance .....                             | 13        |
| 4.3      | Installation .....                                      | 14        |
| 4.3.1    | Balances with draft shield .....                        | 14        |
| 4.3.2    | Balances without draft shield .....                     | 16        |
| 4.4      | Putting into operation .....                            | 17        |
| 4.4.1    | Connecting the balance .....                            | 17        |
| 4.4.2    | Switching on the balance .....                          | 17        |
| 4.4.3    | Leveling the balance .....                              | 18        |
| 4.4.4    | Performing an internal adjustment .....                 | 18        |
| 4.4.5    | Entering / Exiting standby mode .....                   | 18        |
| 4.4.6    | Entering / Exiting power-saving mode .....              | 18        |
| 4.4.7    | Switching off the balance .....                         | 19        |
| 4.5      | Performing a simple weighing .....                      | 19        |
| 4.5.1    | Opening and closing draft shield doors .....            | 19        |
| 4.5.2    | Zeroing the balance .....                               | 19        |
| 4.5.3    | Taring the balance .....                                | 19        |
| 4.5.4    | Performing a weighing .....                             | 19        |
| 4.6      | Transporting, packing, and storing .....                | 20        |
| 4.6.1    | Transporting the balance over short distances .....     | 20        |
| 4.6.2    | Transporting the balance over long distances .....      | 20        |
| 4.6.3    | Packing and storing .....                               | 20        |

|          |                                             |           |
|----------|---------------------------------------------|-----------|
| <b>5</b> | <b>Maintenance</b>                          | <b>21</b> |
| 5.1      | Maintenance tasks .....                     | 21        |
| 5.2      | Cleaning .....                              | 21        |
| 5.2.1    | Disassembling for cleaning.....             | 21        |
| 5.2.1.1  | Balances with draft shield.....             | 21        |
| 5.2.1.2  | Balances without draft shield .....         | 23        |
| 5.2.2    | Cleaning the balance .....                  | 24        |
| 5.2.3    | Putting into operation after cleaning ..... | 24        |
| 5.3      | Service.....                                | 25        |
| <b>6</b> | <b>Technical Data</b>                       | <b>25</b> |
| 6.1      | General data.....                           | 25        |
| <b>7</b> | <b>Disposal</b>                             | <b>26</b> |
| <b>8</b> | <b>Compliance Information</b>               | <b>26</b> |

# 1 Introduction

Thank you for choosing a METTLER TOLEDO balance. The balance combines high performance with ease of use.

## EULA

The software in this product is licensed under the METTLER TOLEDO End User License Agreement (EULA) for Software.

When using this product you agree to the terms of the EULA.

► [www.mt.com/EULA](http://www.mt.com/EULA)

## 1.1 Document purpose

This User Manual provides brief instructions about the first steps to take with the instrument. This ensures safe and efficient handling. Personnel must have carefully read and understood this manual before performing any task.

## 1.2 Further documents and information

This document is available in other languages online.

► [www.mt.com/MR-UM](http://www.mt.com/MR-UM)

Product page:

► [www.mt.com/MR-balances](http://www.mt.com/MR-balances)

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

Search for software:

► [www.mt.com/labweighing-software-download](http://www.mt.com/labweighing-software-download)

Search for documents:

► [www.mt.com/library](http://www.mt.com/library)

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

► [www.mt.com/contact](http://www.mt.com/contact)

## 1.3 Acronyms and abbreviations

| Original term | Explanation                                   |
|---------------|-----------------------------------------------|
| AC            | Alternating Current                           |
| ASTM          | American Society for Testing and Materials    |
| DC            | Direct Current                                |
| EMC           | Electromagnetic Compatibility                 |
| FCC           | Federal Communications Commission             |
| GWP           | Good Weighing Practice                        |
| ID            | Identification                                |
| IP            | Ingress Protection                            |
| LAN           | Local Area Network                            |
| LED           | Light-Emitting Diode                          |
| LPS           | Limited Power Source                          |
| MAC           | Media Access Control                          |
| MT-SICS       | METTLER TOLEDO Standard Interface Command Set |

|      |                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------|
| NA   | Not Applicable                                                                                      |
| OIML | Organisation Internationale de Métrologie Légale<br>(International Organization of Legal Metrology) |
| RM   | Reference Manual                                                                                    |
| SOP  | Standard Operating Procedure                                                                        |
| TDNR | Type Definition Number                                                                              |
| UM   | User Manual                                                                                         |
| USB  | Universal Serial Bus                                                                                |
| USP  | United States Pharmacopeia                                                                          |

## 2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is available online in various languages.
- A printed version of the User Manual is delivered with the instrument.
- The Reference Manual is available online. This manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

### 2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

#### Signal words

|                |                                                                                                                                                         |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DANGER</b>  | A hazardous situation with high risk, resulting in death or severe injury if not avoided.                                                               |
| <b>WARNING</b> | A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.                                                    |
| <b>CAUTION</b> | A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.                                                              |
| <b>NOTICE</b>  | A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data. |

#### Warning symbols

General hazard

Notice

### 2.2 Product-specific safety notes

#### Intended use

This instrument is designed to be used by trained staff. The instrument is intended for weighing purposes.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

#### Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

## Safety notes

### **WARNING**

#### **Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

### **NOTICE**

#### **Damage to the instrument or malfunction due to the use of unsuitable parts**

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

A list of spare parts and accessories can be found in the Reference Manual.

## 3 Design and Function


For further information, consult the Reference Manual (RM).

 [www.mt.com/MR-RM](http://www.mt.com/MR-RM)

### 3.1 Overview analytical balances


0.1 mg

|          |                           |           |                                      |
|----------|---------------------------|-----------|--------------------------------------|
| <b>1</b> | Terminal                  | <b>7</b>  | Top door, draft shield               |
| <b>2</b> | Drip tray                 | <b>8</b>  | Side door, draft shield (right/left) |
| <b>3</b> | Weighing pan              | <b>9</b>  | Handle, side door                    |
| <b>4</b> | Front panel, draft shield | <b>10</b> | Back panel, draft shield             |
| <b>5</b> | Handle, top door          | <b>11</b> | Leveling feet                        |
| <b>6</b> | QuickLock, doors/panel    |           |                                      |

### 3.2 Overview precision balances, with draft shield


1 mg

|          |                           |           |                                      |
|----------|---------------------------|-----------|--------------------------------------|
| <b>1</b> | Terminal                  | <b>7</b>  | Top door, draft shield               |
| <b>2</b> | Drip tray                 | <b>8</b>  | Side door, draft shield (right/left) |
| <b>3</b> | Weighing pan              | <b>9</b>  | Handle, side door                    |
| <b>4</b> | Front panel, draft shield | <b>10</b> | Back panel, draft shield             |
| <b>5</b> | Handle, top door          | <b>11</b> | Leveling feet                        |
| <b>6</b> | QuickLock, doors/panel    |           |                                      |


### 3.3 Overview precision balances, without draft shield

0.01 g


|   |                          |   |               |
|---|--------------------------|---|---------------|
| 1 | Terminal                 | 4 | Weighing pan  |
| 2 | Draft-protection element | 5 | Leveling feet |
| 3 | Weighing pan support     |   |               |

### 3.4 Overview terminal

| Name                                                                                                                 | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <b>Standby / Power-saving mode</b> | By tapping  , the balance goes into standby mode. By tapping and holding  , the balance goes into power-saving mode. To switch the balance completely off, it must be unplugged from the power supply.<br><br><b>Note</b><br>Do not disconnect the balance from the power supply unless the balance is not used for an extended period of time. After switching on the instrument, it must warm up before giving accurate results. |
|  Tare                               | Tares the balance.<br><br>This function is used when the weighing process involves containers. After taring the balance, the screen shows <b>Net</b> which indicates that all displayed weights are net.                                                                                                                                                                                                                                                                                                                                                                                              |
|  Zero                               | Zeroes the balance.<br><br>The balance must always be zeroed before starting the weighing process. After zeroing, the balance sets a new zero point.                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|  Home                               | To return from any menu level to the main weighing screen.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

### 3.5 Overview interface connections


|          |                           |          |                          |
|----------|---------------------------|----------|--------------------------|
| <b>1</b> | RS232C serial interface   | <b>4</b> | USB-A port               |
| <b>2</b> | Slot for anti-theft cable | <b>5</b> | USB-C port               |
| <b>3</b> | Ethernet port (LAN)       | <b>6</b> | Socket for AC/DC adapter |

### 3.6 Components description


#### 3.6.1 Draft shield

The draft shield protects the weighing area against environmental impacts like drafts or moisture. The side doors and the top door can be opened manually.


### 3.6.2 Weighing pan

The weighing pan is the load receptor that serves to accommodate the weighing item.


### 3.6.3 Drip tray

The drip tray is positioned below the weighing pan. The primary purpose of the drip tray is to ensure quick cleaning of the balance.


### 3.6.4 Door handle

The door handles are mounted on the draft shield doors. The handles are used to manually open the side doors and the top door of the draft shield.


### 3.6.5 Leveling feet

The balance stands on height-adjustable feet. These feet are used to level the balance.


**3.6.6 Terminal**

The balance terminal has a 4.3-inch touch-sensitive display. The terminal and the platform are protected by a replaceable cover.


**3.6.7 QuickLock for draft shield**

The QuickLock for the draft shield is used to secure the draft shield to the platform.


**3.6.8 QuickLock for doors and front panel**

Depending on the position, the QuickLock is used to lock/unlock the top door, the side doors, and the front panel of the draft shield.

**3.6.9 Release button for back panel**


The release button is used to lock/unlock the back panel of the draft shield.

## 3.7 User interface

### 3.7.1 Main sections at a glance


The main weighing screen (1) is the central navigation point where all the menus and settings can be found. The sections **Balance menu** (2), **Main configuration** (3), and the applications section (4) open when tapping the corresponding icon or tab.

#### See also

Main weighing screen ► Page 12

3.7.2 Main weighing screen

|    | Name                      | Description                                                                                                                                                                          |
|----|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Weighing results          | Shows the results of the current weighing process.                                                                                                                                   |
| 2  | Level indicator           | Indicates if the balance is leveled (green) or not (red).                                                                                                                            |
| 3  | Applications              | Accesses available applications: <b>Weighing, Adjustments, Tests.</b>                                                                                                                |
| 4  | Readability and capacity  | Shows the readability and the capacity of the balance.                                                                                                                               |
| 5  | Additional information    | Shows additional information about the current activity.<br>Example: current weighing value in another unit                                                                          |
| 6  | Information and warnings  | Shows current information, warnings, and error messages.                                                                                                                             |
| 7  | Functions area            | Shows the active functions according to the settings of the current weighing application.                                                                                            |
| 8  | <b>Publish</b> button     | Publishes the results according to the settings of the current weighing application.<br><br>Depending on the selected weighing application, the button can have different functions. |
| 9  | Action bar                | Contains actions referring to the current weighing application.                                                                                                                      |
| 10 | <b>Balance menu</b>       | Accesses the balance properties.                                                                                                                                                     |
| 11 | SmartTrac                 | Used as a weighing aid to define a target weight with upper and lower tolerances.                                                                                                    |
| 12 | <b>Main configuration</b> | Accesses the configuration options for the current weighing application.                                                                                                             |

## 4 Installation and Putting into Operation

### 4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

#### Requirements of the location

Place indoors on stable table



Ensure sufficient spacing

Level the instrument

Provide adequate lighting


Avoid direct sunlight


Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Take into account the environmental conditions. See "Technical Data".

Sufficient spacing for balances: > 15 cm all around the instrument


### 4.2 Unpacking the balance

Check the package, the packaging elements and the delivered components for damages. If any components are damaged, please contact your METTLER TOLEDO service representative.


#### Note

Depending on the balance model, the packaging elements and the components may look different.


- 1 Open the box and lift the package out using the lifting strap (1).
- 2 Place the package on a level surface with the inscription BOTTOM (2) facing downwards.


- 3 Open the lifting strap (1) and remove the User Manual (3).

- 4 Remove the upper part of the package and unpack the drip tray (4).

- 5 Carefully unpack the balance (5) and all other items.  
6 Remove the protective bag.  
7 Keep the protective cover installed on the platform and on the terminal.  
8 Store all parts of the packaging in a safe place for future use.  
⇒ The balance is ready for installation.

## 4.3 Installation

### Note

Depending on the balance model, the components may look different.

### 4.3.1 Balances with draft shield

### CAUTION


#### Injury due to sharp objects or broken glass

- Instrument components, e.g., glass, can break and lead to injuries.
- Always proceed with focus and care.


### Note

Skip steps 1 – 4 when unpacking the balance for the first time, or if the draft shield is already assembled.


- 1 Assemble the draft shield: Turn the QuickLock (1, right, left) and slide in the top door (2).

- 2 Slide in the side door (3) (right, left).


- 3 Attach the front panel (4), then turn the QuickLock (1, right, left) to hold the panel in place.

- 4 Attach the back panel (5). Make sure the release button (6) snaps in.


➡ The draft shield is assembled.

- 5 Secure the draft shield (7) to the platform (9) by turning the QuickLock (8).

**Note**


To protect your balance, keep the protective cover installed on the platform (9).

- 6 Insert the drip tray (10).


**Note**

The ring-shaped elevation (11) must point upwards. It serves as draft protection.

- 7 Install the weighing pan (12).

➡ The balance is ready for use.

#### 4.3.2 Balances without draft shield


- 1 Place the draft-protection element (1) on top of the platform (2).

**Note**

To protect your balance, keep the protective cover installed on the platform (2).

- 2 Place the weighing pan support (3) on top of the platform (2).
- 3 Place the weighing pan (4) on the weighing pan support (3).

➡ The balance is ready for use.

## 4.4 Putting into operation

### 4.4.1 Connecting the balance


#### **WARNING**

##### **Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
  - 2 Insert the plug of the AC/DC adapter (1) into the power socket of the instrument.
  - 3 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- ⇒ The balance automatically switches on.

#### **Note**

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

#### **See also**

🔗 General data ▶ Page 25

### 4.4.2 Switching on the balance

When connected to the power supply, the balance automatically switches on.

#### **EULA (End User License Agreement)**

When the balance is switched on the first time, the EULA (End User License Agreement) appears on the screen.

- 1 Read the conditions.
  - 2 Tap **I accept the terms in the license agreement** and confirm with **✓ OK**.
- ⇒ The main weighing screen appears.

#### **Acclimatization and warm up**

Before the balance gives reliable results, it must:

- acclimatize to the room temperature
- warm up by being connected to the power supply

The acclimatization time and warm-up time for balances are available in "General data".

#### **Note**

When the balance is exiting standby, it is ready immediately.

### See also

- 🔗 General data ▶ Page 25
- 🔗 Entering / Exiting standby mode ▶ Page 18
- 🔗 Entering / Exiting power-saving mode ▶ Page 18
- 🔗 Switching off the balance ▶ Page 19

#### 4.4.3 Leveling the balance


Exact horizontal and stable positioning are essential for repeatable and accurate weighing results.

If the balance is out of level, the level indicator on the main screen turns red.

- 1 On the main weighing screen, tap .
  - ➔ The dialog **Leveling aid** opens.
- 2 Turn both leveling feet (1) as instructed on the display until the dot is in the center of the level indicator.

Alternative access to the dialog **Leveling aid**:

≡ **Navigation:** ≡ **Balance menu** >  **Leveling aid**

#### 4.4.4 Performing an internal adjustment

≡ **Navigation:** ▼ >  **Applications** >  **Adjustments**

- **Adjustments** is set to **Internal**.

- 1 Option 1: On the main weighing screen, tap  **Adjust**.  
Option 2: Open the applications section, tap  **Adjustments**, select the adjustment, and tap ▶ **Start**.
  - ➔ The adjustment is executed.
  - ➔ The adjustment results appear.
- 2 Tap  **Finish**.
  - ➔ The balance is ready.

#### 4.4.5 Entering / Exiting standby mode

- 1 To enter standby mode, short press .
  - ➔ The display is blue. A QR code for more information about the balance is shown.
- 2 To exit standby mode, short press .
  - ➔ The display is turned on.

#### 4.4.6 Entering / Exiting power-saving mode

- 1 To enter power-saving mode, long press  (more than two seconds).
  - ➔ The display is dark. The balance is in power-saving mode.
- 2 To exit power-saving mode, long press .
  - ➔ The balance is switched on.

#### **Note**

We recommend configuring power-saving times. When the balance exits power-saving mode automatically at the defined time, the balance is ready for use immediately.

If the power-saving mode is terminated manually, the balance must warm up before it can be used.

### See also

- 🔗 General data ▶ Page 25

#### 4.4.7 Switching off the balance

To completely switch off the balance, it must be disconnected from the power supply. By pressing , the balance goes only into standby mode or into power-saving mode.

##### Note

When the balance has been completely switched off for some time, it must warm up before it can be used.

##### See also

-  Switching on the balance ► Page 17
-  Entering / Exiting standby mode ► Page 18
-  Entering / Exiting power-saving mode ► Page 18


### 4.5 Performing a simple weighing

##### Note

A balance with draft shield is used to explain the procedure. For balances without a draft shield, skip the instructions steps concerning the draft shield.

#### 4.5.1 Opening and closing draft shield doors

- Open the door manually with the door handle (1).


#### 4.5.2 Zeroing the balance

- 1 Open the draft shield.
- 2 Clear the weighing pan.
- 3 Close the draft shield.
- 4 Press  to zero the balance.
  - ⇒ The balance is zeroed.

#### 4.5.3 Taring the balance

If a sample vessel is used, the balance must be tared.

- The balance is zeroed.
- 1 Place the sample vessel (1) on the weighing pan (2).
- 2 Press  to tare the balance.
  - ⇒ The balance is tared. The icon  appears.

#### 4.5.4 Performing a weighing

- 1 Open the draft shield.
- 2 Place the weighing object into the sample vessel.

- 3 Close the draft shield.  
➡ The result is displayed.
- 4 Optional, if a printer is connected: Tap  to print the weighing result.


## 4.6 Transporting, packing, and storing

### 4.6.1 Transporting the balance over short distances

- 1 Disconnect the AC/DC adapter and unplug all interface cables.
- 2 Hold the balance with both hands and carry it in horizontal position to the target location. Consider the requirements of the location.

If you want to put the balance into operation, proceed as follows:

- 1 Connect in reverse order.
- 2 Give the balance sufficient time to warm up.
- 3 Level the balance.
- 4 Perform an internal adjustment.

#### See also

-  Selecting the location ▶ Page 13
-  Switching on the balance ▶ Page 17
-  Leveling the balance ▶ Page 18
-  Performing an internal adjustment ▶ Page 18

### 4.6.2 Transporting the balance over long distances

METTLER TOLEDO recommends using the original packaging for transportation or shipment of the balance or balance components over long distances. The elements of the original packaging are developed specifically for the balance and its components and ensure maximum protection during transportation.

#### See also

-  Unpacking the balance ▶ Page 13

### 4.6.3 Packing and storing

#### Packing the balance

Store all parts of packaging in a safe place. The elements of the original packaging are developed specifically for the balance and its components, and ensures maximum protection during transportation and storage.

#### Storing the balance

Only store the balance under the following conditions:

- Indoor and in the original packaging
- According to the environmental conditions, see "Technical Data"

#### Note

When storing for longer than two weeks, the battery (capacitor) may become empty (only date and time get lost).

#### See also

-  Technical Data ▶ Page 25

## 5 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

For further information, consult the Reference Manual (RM).

► [www.mt.com/MR-RM](http://www.mt.com/MR-RM)

### 5.1 Maintenance tasks

| Maintenance action                                                                                                                                       | Recommended interval                                                                                                                                                                          | Remarks                                       |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| Performing an internal adjustment                                                                                                                        | <ul style="list-style-type: none"> <li>Daily</li> <li>After cleaning</li> <li>After leveling</li> <li>After changing the location</li> </ul>                                                  | see "Performing an internal adjustment"       |
| Performing routine tests (eccentricity test, repeatability test, sensitivity test).<br>METTLER TOLEDO recommends to at least perform a sensitivity test. | <ul style="list-style-type: none"> <li>After cleaning</li> <li>After assembling the balance</li> <li>After a software update</li> <li>Depending on your internal regulations (SOP)</li> </ul> | see "Tests" in the Reference Manual           |
| Cleaning                                                                                                                                                 | <ul style="list-style-type: none"> <li>After every use</li> <li>Depending on the degree of pollution</li> <li>Depending on your internal regulations (SOP)</li> </ul>                         | see "Cleaning"                                |
| Updating the software                                                                                                                                    | <ul style="list-style-type: none"> <li>Depending on your internal regulations (SOP).</li> <li>After a new software release.</li> </ul>                                                        | see "Software update" in the Reference Manual |

#### See also

🔗 Performing an internal adjustment ► Page 18

🔗 Cleaning ► Page 21

## 5.2 Cleaning

### 5.2.1 Disassembling for cleaning

#### Note

Depending on the balance model, the components may look different.

#### Note

In most cases, it is not necessary to remove the protective cover to clean the balance.

#### 5.2.1.1 Balances with draft shield


#### CAUTION

##### Injury due to sharp objects or broken glass


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Remove the weighing pan (1) and the drip tray (2).


- 2 Turn the QuickLock (3, right, left) and pull the side door (4) towards the front to remove it (right, left).

- 3 Tilt the front panel (5) towards the front and lift it upwards to remove it.


- 4 Pull the top door (6) towards the front to remove it.


- 5 Press the release button (7) and tilt the back panel (8) to remove it.

**Note**


Optional, if required: Remove the protective cover for cleaning as described below.

- 6 Open the QuickLock (9) and remove the draft shield (10).

- 7 Remove the protective cover (11) from the platform (12).


### 5.2.1.2 Balances without draft shield


- 1 Remove the draft-protection element (1).
- 2 Remove the weighing pan (2).
- 3 Remove the weighing pan support (3).

**Note**

Optional, if required: Remove the protective cover for cleaning as described below.

- 4 Remove the screws (4) to remove the EMC plate (5).
- 5 Remove the protective cover (6) from the platform (7).

### 5.2.2 Cleaning the balance

#### NOTICE

##### **Damage to the instrument due to inappropriate cleaning methods**

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

#### **Cleaning around the balance**

- Remove any dirt or dust around the balance and avoid further contaminations.

#### **Cleaning the terminal**

- Clean the terminal with a damp cloth or a tissue and a mild cleaning agent.

#### **Cleaning the removable parts**

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent or clean in a dishwasher up to 80 °C.

#### **Cleaning the balance**

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the balance.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

### 5.2.3 Putting into operation after cleaning

- 1 Reassemble the balance.
- 2 Check that the draft shield doors (top, sides) open and close normally (if applicable).
- 3 Reconnect the balance to the AC/DC adapter.
- 4 Check the level status, level the balance if necessary.

- 5 Respect the warm-up time specified in the "Technical Data".
- 6 Perform an internal adjustment.
- 7 Perform a routine test according to the internal regulations of your company. METTLER TOLEDO recommends performing a sensitivity test after cleaning the balance.
- 8 Press **→0←** to zero the balance.
  - ⇒ The balance is ready for use.

#### See also

- 🔗 Leveling the balance ► Page 18
- 🔗 Technical Data ► Page 25
- 🔗 Performing an internal adjustment ► Page 18

## 5.3 Service

Regular servicing by an authorized service technician ensures reliability for years to come. Contact your METTLER TOLEDO representative for details about the available service options.

## 6 Technical Data

### 6.1 General data

#### Power supply

|                            |                                                                             |
|----------------------------|-----------------------------------------------------------------------------|
| AC/DC adapter:             | Input: 100 – 240 V AC ± 10%, 50 – 60 Hz, 0.5 A<br>Output: 12 V DC, 1 A, LPS |
| Balance power consumption: | 12 V DC, 0.8 A                                                              |
| Polarity:                  | ⬡—●—⬡                                                                       |

#### Protection and standards

|                          |                                                      |
|--------------------------|------------------------------------------------------|
| Overvoltage category:    | II                                                   |
| Degree of pollution:     | 2                                                    |
| Ingress protection code: | IP43 (balances with readability of 0.01 g or higher) |

#### Note

Stated IP is only achieved when the balance is ready for operation. The protective cover must be installed, and the caps must cover the interface connections.

|                               |                                   |
|-------------------------------|-----------------------------------|
| Standards for safety and EMC: | See Declaration of Conformity     |
| Range of application:         | Use only indoors in dry locations |

#### Environmental conditions

The limit values apply when the balance is used under the following environmental conditions:

|                              |                                                                                                                                                                                                                      |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Height above mean sea level: | Up to 5000 m                                                                                                                                                                                                         |
| Ambient temperature:         | +10 – +30 °C                                                                                                                                                                                                         |
| Temperature change, max.:    | 5 °C/h                                                                                                                                                                                                               |
| Relative humidity:           | 30 – 70%, non-condensing                                                                                                                                                                                             |
| Acclimatization time:        | Recommendation: Up to <b>4 hours</b> for precision balances, or up to <b>8 hours</b> for analytical balances. These values apply after placing the balance in the same location where it will be put into operation. |

#### Note

The acclimatization time depends on the readability of the balance, and on the environmental conditions.

Warm-up time:

At least **30 minutes** for precision balances, or **60 minutes** for analytical balances. These values apply after connecting the balance to the power supply, or after exiting power-saving mode. When switched on from standby, the balance is ready for operation immediately.

The balance can be used under the following environmental conditions. However, the weighing performances of the balance may be outside the limit values:

Ambient temperature: +5 °C – +40 °C

Relative humidity: 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

The balance can be disconnected and stored in its packaging under the following conditions:

Ambient temperature: -25 – +70 °C

Relative humidity: 10 – 90%, non-condensing

## 7 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

## 8 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► [www.mt.com/ComplianceSearch](http://www.mt.com/ComplianceSearch)

For further information, consult the Reference Manual (RM).

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► [www.mt.com/contact](http://www.mt.com/contact)

## 1.3 Akronymmer og forkortelser

| Originalt udtryk | Oversat udtryk | Forklaring                                 |
|------------------|----------------|--------------------------------------------|
| AC               |                | Alternating Current                        |
| ASTM             |                | American Society for Testing and Materials |
| DC               |                | Direct Current                             |
| EMC              |                | Electromagnetic Compatibility              |
| FCC              |                | Federal Communications Commission          |
| GWP              |                | Good Weighing Practice                     |
| ID               |                | Identification                             |
| IP               |                | Ingress Protection                         |
| LAN              |                | Local Area Network                         |
| LED              |                | Light-Emitting Diode                       |
| LPS              |                | Limited Power Source                       |
| MAC              |                | Media Access Control                       |

|         |                                                                                                  |
|---------|--------------------------------------------------------------------------------------------------|
| MT-SICS | METTLER TOLEDO Standard Interface Command Set                                                    |
| NA      | Not Applicable                                                                                   |
| OIML    | Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology) |
| RM      | Reference Manual                                                                                 |
| SOP     | Standard Operating Procedure                                                                     |
| TDNR    | Type Definition Number                                                                           |
| UM      | User Manual<br>(Brugervejledning)                                                                |
| USB     | Universal Serial Bus                                                                             |
| USP     | United States Pharmacopeia                                                                       |

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© 07/2024 METTLER TOLEDO. All rights reserved.