English  
中文  
Deutsch  
Français  
Español  
Русский

**METTLER TOLEDO**

**DX205-SC**

**Reference Sensor**

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**Referenzsensor**

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Standard reference sensor for ion-selective measurement with movable PTFE-sleeve diaphragm

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## DX205-SC

### Overview

- |   |   |
|---|---|
| 1 Sleeve diaphragm (Movable PTFE sleeve)                                | 2 Ag/AgCl reference system                    |
| 3 Bridge electrolyte chamber, 0.5 mol/L Diisopropylamine - ISA solution | 4 Shaft (glass)                               |
| 5 Safe lock   | 6 Filling aperture for the bridge electrolyte |
| 7 Electrode head (S7 connector)   |   |

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## 1 Safety notes

### Note

#### Electrolyte causes skin irritation!

Diisopropylamin causes skin irritation. For more information on the handling with the electrolyte, refer to the related safety data sheet.

### Note

#### Damage to the Sensor!

The sensor will be damaged if it is allowed to dry out.

Make sure that the sensor is kept in a suitable liquid at all times.

Suitable eye protection such as goggles should be worn.

Use appropriate gloves when handling chemicals or hazardous substances, checking their integrity before use.

## 2 Characteristics

- Temperature range: 0 - 60 °C (32-132°F).
- Movable PTFE reference junction.

## 3 Start-up maintenance

### Note

- White deposits in the sensor have no influence on the measurement result!
- The electrode is supplied without electrolyte filling.

Fill up the sensor with the recommended reference electrolyte as indicated below:

- 1 Remove the nipple from the filling aperture (6).
- 2 Fill up with reference electrolyte (0.5 mol/L DIPA - HCl - ISA) to approx. 10 mm below the filling aperture (6).
- 3 Check by hand that the sleeve diaphragm (1) is smooth running.
- 4 Rinse the sensor with distilled water.

**Note:** If the electrode has been stored for a long time without bridge electrolyte (especially when using the electrode for the first time), the freshly refilled electrode must be conditioned for at least 24h in the reference electrolyte.

⇒ The sensor is ready for use.

## 4 Storage

- 1 Rinse and clean the sensor after use and place the wetting cap over the sleeve diaphragm (1).
- 2 If necessary, refill the sensor with the recommended reference electrolyte (min level = 1/3).
- 3 Make sure the filling aperture (6) is closed.
- 4 If the sensor remains in the analysis beaker, make sure the sleeve diaphragm (1) is fully immersed in liquid (e.g. reference electrolyte in the analysis beaker).

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© Mettler-Toledo AG 09/2013
**Mettler-Toledo AG, Analytical**