METTLER TOLEDO

InPro 7250 Series

**Inductive Conductivity Sensors
Induktive Leitfähigkeitsmesszellen
Sondes de conductivité inductive**

**Instruction Manual
Bedienungsanleitung
Instructions d'utilisation**

English **Page 3**

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InPro 7250 Series

Instruction Manual

© It is forbidden to reprint this Instruction Manual in whole or part. No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other recordings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH - 8902 Urdorf, Switzerland.

All rights reserved, in particular reproduction, translation and patenting/registration.

Contents

1	Introduction.....	5
2	Important notes	6
2.1	Notes on operating instructions.....	6
2.2	Intended use.....	6
2.3	Safety instructions	7
2.4	Correct disposal of the sensor.....	9
2.5	Examples of some typical applications	9
2.6	Use in Ex zones	9
2.7	Ex classification ATEX	10
2.7.1	Introduction	10
2.7.2	Rated data	10
2.7.3	Special conditions.....	11
2.8	Ex classification FM approved	12
2.8.1	Introduction	12
2.8.2	Rated data	12
2.8.3	Special conditions.....	12
2.8.4	Applied standards	13
3	Product description.....	14
3.1	General information.....	14
3.2	Principle	14
3.3	Sensor identification	14
3.4	Sensor applications.....	15
4	Installation, operation and maintenance	16
4.1	Sensor installation.....	16
4.2.	Connection of the sensor.....	18
4.3	Flange installation.....	18
4.4	Bushing installation.....	19
4.5	Immersion installation with an InDip 550 Ind	19
4.6	Sanitary installation.....	20
4.7	Environmental protection.....	20
5	Product specifications.....	21
5.1.	Sensor specifications.....	21
6	Ordering information	22
6.1.	Sensors	22
6.2	Spare parts.....	22
6.3	Process Connections and Accessories	22
6.4	Transmitter M700(x), module Cond Ind 7700(x).....	23
6.5	Transmitter M400, (4-wire Transmitter).....	23
6.6	Transmitter M400, (2-wire Transmitter).....	23
	Appendix	67
	Sensor permissible pressure and temperature diagram	67

InPro and InDip are trademarks
of the METTLER TOLEDO Group.

Viton is a trademark
of DuPont Performance Elastomers LLC.

All other trademarks are the property
of their respective holders.

1 Introduction

Thank you for buying the **InPro 7250™ Series sensor from METTLER TOLEDO.**

The construction of the InPro 7250 Series sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment. **Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Important notes

2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 7250 Series sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to safety instructions and warnings of potential danger which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies additional information and instructions which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended use

The InPro 7250 Sensor Series is intended to be used to measure conductivity and concentration of liquids, in accordance with the data in this instruction manual. Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with METTLER TOLEDO Process Analytics in the event of any uncertainties.

2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The conductivity sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see ordering information, "Section 6").
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damage caused by unauthorized modifications. The risk is borne entirely by the user.
- When using cables for sensor connection not supplied and recommended by the manufacturer the following standards need to be considered:
 - US National Electrical Code® (ANSI/NFPA 70 [NEC®]), where applicable.
 - Canadian Electrical (CE) Code® (CEC Part 1, CAN/CSA-C22.1), where applicable.

-
- For guidance on US installations, see ANSI/ISA-RP12.06.01, Installation of Intrinsically Safe Systems for Hazardous (Classified) Locations.
 - Control room equipment connected to intrinsically safe associated apparatus in the US should not use or generate more than the specified Um of the associated apparatus.
 - Care must be taken during installation to avoid impacts or friction that could create an ignition source.
 - Tampering and replacement with non-factory components may adversely affect the safe use of the system.
 - Insertion or withdrawal of removable electrical connectors or modules is to be accomplished only when the area is known to be free of flammable vapors.
 - WARNING – INTRINSICALLY SAFE APPARATUS CAN BE A SOURCE OF IGNITION IF INTERNAL SPACINGS ARE SHORTED OR CONNECTIONS OPENED.
 - WARNING – DO NOT OPEN WHEN AN EXPLOSIVE ATMOSPHERE IS PRESENT.
 - WARNING – SUBSTITUTION OF COMPONENTS MAY IMPAIR INTRINSIC SAFETY.
 - WARNING – SUBSTITUTION OF COMPONENTS MAY IMPAIR SUITABILITY OF THE EQUIPMENT
 - WARNING – FOR CONNECTION ONLY TO NON-FLAMMABLE PROCESSES.
 - WARNING – TO MAINTAIN THE ENCLOSURE IP66 RATING, THE CONNECTOR MUST BE FULLY ENGAGED.

2.4 Correct disposal of the sensor

When the sensor is finally removed from service, observe all local environmental regulations for proper disposal. (see also 4.7 Environmental protection)

2.5 Examples of some typical applications

Below is a list of examples of typical fields of application for the conductivity sensors. This list is not exhaustive.

Measurement in liquids:

- Industrial waste water
- Salinity and brine measurements
- Steel pickling
- Scrubbing towers
- Ion exchange regeneration
- Plating baths
- Rinse water
- Caustic metal cleaning
- Textile measurements
- Oleum measurements
- CIP skids
- Acid and base concentration

2.6 Use in Ex zones

Attention!

For an installation in Ex zones please read the guidelines following hereafter:

Ex classification ATEX:

 II 1/2G Ex ia IIC T6...T3 Ga/Gb

Number of the test certificate:

SEV 15 ATEX 0123X
IECEX SEV 15.0011X
CML 22 UKEX 2186X

Ex classification FM approved:

 IS/I,II,III/1/ABCDEFG/T6 Ta = 60 °C

Designation and number of the declaration:

Project ID: 3021227

FM certificate number: FM16US0034X

FM18CA0021X

2.7 Ex classification ATEX

2.7.1 Introduction

According to RL2014/34/EU¹⁾ (ATEX 114) Appendix I, InPro 725X/**/* conductivity sensors are devices group II, category 1/2G and according to RL 99/92/EG²⁾ (ATEX 137) may be used in zones 0/1 or 1/2 or 1 or 2 and gas groups IIA, IIB and IIC that are potentially explosive due to combustible substances in above stated temperature classes.

For use / installation, the requirements of EN 60079-14 must be observed.

According to RL2014/34/EU¹⁾ (ATEX 114) Appendix I, InPro 725X/**/* conductivity sensors are devices group II, category 1/2D and according to RL 99/92/EG²⁾ (ATEX 137) may also be used in zones 20/21 or 21 or 22 that contain combustible dusts.

The conductivity measuring circuit and temperature measuring circuit are part of a common intrinsically safe system and is for operation connected to a separately certified transmitter.

The conductivity measuring circuit and temperature measuring circuit as part of an intrinsically safe system is isolated from the not-intrinsically safe electric circuits up to a maximum rated voltage of 375 V and from grounded parts up to a maximum rated voltage of 30 V.

The sensors can be installed with different process adapters (flange, threaded bushes) permanently in pipes or tanks.

2.7.2 Rated data

Conductivity measuring circuit, and temperature measuring circuit 1:

In type of protection intrinsic safe Ex ia IIC, only for connection to a certified intrinsically safe circuit.

Maximum values:

$U_i \leq 16 \text{ V}$, $I_i \leq 150 \text{ mA}$, $P_i \leq 155 \text{ mW}$

$L_i = 0$ (the internal inductance is ineffective towards the outside)

$C_i = 900 \text{ pF}$ (effective internal capacitance)

Note:

The above maximum values are each the total of all individual circuits of the associated intrinsically safe power supply and transmitter.

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR).

2.7.3 Special conditions

1. The relationship between the maximum permissible ambient or media temperature and temperature class is shown in the following table :

Temperature class	Max. environment resp. media temperature
T6	68 °C
T5	80 °C
T4	108 °C
T3	130 °C


2. The capacitance and inductance of the connecting cable has to be considered.
3. The InPro 725X/**/* conductivity sensors may only be used in suitable process terminals of METTLER TOLEDO or other manufacturers in potentially explosive atmospheres.
4. The independent process terminal used for installation of the conductivity sensors must be connected to the equipotential bonding system of the installation.
5. The independent process terminal used for installation of the conductivity sensors must be included in the recurring pressure test of the installation if necessary.
6. The minimum conductivity of the media for safe working in potentially explosive atmospheres must be higher than 1 nS/cm.

2.8 Ex classification FM approved

2.8.1 Introduction

The following FM control drawing and the standards listed in section 2.8.4 must be observed, where applicable:

2.8.2 Rated data

For rated data please observe section 2.7.2

2.8.3 Special conditions

In type of protection intrinsically safe apparatus, the analytical probes (InPro 7000 Series) equipment is designated with the following specific conditions of use.

1. Process temperature no greater than +130 °C.
2. Maximum permissible working pressure is 15 barg (218 psig).
3. Potential Electrostatic Charging Hazard – To prevent the risk of electrostatic sparking, the non-metallic surface should only be cleaned with a damp cloth.

2.8.4 Applied standards

United States Standards

Title	Number	Issue Date
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

Canadian Standards

Title	Number	Issue Date
Explosive Atmospheres – Part 0: Equipment – General Requirements	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	CAN/CSA-C22.2 No. 61010-1	R2017

3 Product description

3.1 General information

The conductivity sensors InPro 7250 Series with integrated temperature sensor are intended to be used to measure conductivity and concentrations of liquids. The sensor is designed for the measurement of medium and high conductivity values.

3.2 Principle

The sensors InPro 7250 Series measure an induced current in a loop of solution.

Two toroidally wound coils are encapsulated in close proximity within the sensor which is immersed in the solution. An AC signal, applied to one toroidal coil, induces a current in the second coil which is directly proportional to the conductance of the solution.

3.3 Sensor identification

Type designation

A sensor can be identified by the article number and the model name on the label located on the cable.

InPro 7250 xx/yy/zz

xx: Type of sensor

ST = Sensor for temperature up to 100 °C (212 °F), without Ex approval

HT = Sensor for temperature up to 180 °C, (356 °F) with Ex approval

PFA = Sensor in PFA* for Temperatures up to 125°C (257°F)

yy: Type of the integrated temperature sensor

zz: Cable length

Serial number

The serial number is located on the end of the mounting thread.

Sensor certificate number

The labels located on the sensor cable provide the information about the sensor certifications.

*PFA: Perfluoralkoxy-Copolymer


Figure 1: Sensor identification

3.4 Sensor applications

InPro 7250 ST

This sensor may be used for most routine applications like conductivity measurement in industrial wastewater systems. For applications at higher temperatures and for use in hazardous areas the sensor InPro 7250 HT is the preferred alternative.

InPro 7250 HT

The InPro 7250 Series sensor is suitable for the majority of all electrodeless conductivity applications. Its compact size enables it to be mounted in a multitude of methods, including insertion (flange, bushing), and immersion. Typical applications include salinity and brine measurements, steel pickling, scrubbing towers, ion exchange regeneration, plating baths, rinse water, caustic metal cleaning, and textile measurements in scouring, mercerizing, and carbonizing baths.

InPro 7250 PFA

This sensor is particularly suited for heavy-duty chemical applications such as Oleum (highly concentrated sulfuric acid), as well as applications in the pharmaceutical and biopharmaceutical industry (CIP skids, acid and base concentration).

Notes:

- PEEK is a thermoplastic material with excellent strength and chemical resistance properties over a wide range of process temperatures and pressures.
- PEEK material displays excellent chemical resistance to most aqueous solutions of acids, bases, and salts. It is also excellent for organic solvents such as toluene, ethyl acetate, acetone, gasoline, and carbon tetrachloride. It is not recommended for sulfuric or nitric acid solutions above 70 %, nor is it recommended for Oleum applications.
- PFA is a type of fluorocarbons with similar properties as PTFE. PFA is resistant against nearly all chemicals and has a high temperature resistance. The PFA quality used here does not contain glass fibres (virgin PFA).

4 Installation, operation and maintenance

4.1 Sensor installation

The sensor shaft is equipped with a conical G 3/4" mounting thread. The installation of the sensor is made by threaded bushes or flanges. The process area is sealed perfectly with flat gasket and an O-ring.

Note: When mounted in-line, the sensor should be installed centered to avoid pipe wall effects. If there is less than 30 mm (1.18 inch) spacing between sensor and pipe, a product calibration with a sensor built in the process is recommended to meet the specified accuracies.


Figure 2: Sensor overview

Note: Proper installation of the sensor is important for efficient and accurate operation. For all applications and sensor configurations, mounting arrangements must be located so that:







- Sample at the sensing area is representative of the solution.
- Solution circulates actively and continuously past the sensing area.
- Position and orientation of the sensor does not trap air bubbles within sensing area.
- Deposits of sediment or other foreign material do not accumulate within the sensing area. If cable is installed in a conduit (recommended), either flexible conduit should be used or some other provision made for removal of sensor from the process.

**Notes:**

- The sensors of the InPro 7250 Series are marked with an arrow, which indicates the flat (upstream) face of the sensor, for correct installation in a pipe or a tank. The marking is at the upper end of the sensor mounting thread (see **Figure 2**).

**Warning!**

Incorrect O-ring position or a damaged O-ring could result in process leakage, causing injury to personnel.

Recommended**Not recommended!****Air entrapment****Solid deposits****Figure 3:**

Installation examples of the InPro7250 sensor

4.2. Connection of the sensor

Cable configuration

	InPro 7250 ST / PFA	InPro 7250 HT
receive hi	black	red
receive lo	red	yellow
send lo	brown	violet
send hi	blue	black
RTD	green	green
RTD	white	white
RTD sense	grey	grey
shield	green/yellow	green/yellow

Transmitter connection

The sensors InPro 7250 Series are suitable for connection to the following METTLER TOLEDO transmitters:

M 700(X), Module Cond Ind 7700(X):

Transmitter, modular measuring system

M400, Type 1 Cond Ind:

Transmitter, 4-wire

M400 2XH Cond Ind:

Transmitter, 2-wire

Note: In order to connect the sensor with the transmitter, please consider the instruction manual of the METTLER TOLEDO transmitter.

4.3 Flange installation

METTLER TOLEDO flanges are used for permanent installations in pipelines and tanks. Flanges are used with electrodeless conductivity sensors in systems using DN50 or 2-inch and larger process piping. Flanges are provided complete with locknut.


Figure 4: Flange installation

4.4 Bushing installation

Supplied METTLER TOLEDO bushings are used for installations in pipelines and tanks and are provided complete with locknut.


Figure 5: Bushing installation

4.5 Immersion installation with an InDip 550 Ind

For open tanks or inlet channel the sensor can be installed in an immersion housing InDip™ 550 Ind.


Figure 6: Installation with an immersion housing «InDip 550Ind»

4.6 Sanitary installation

METTLER TOLEDO offers a dairy adapter according to DIN 11851 (DN50) as well as an aseptic adapter DIN 11864-1 (DN50) for the inductive sensors of the InPro 7250 series.


Figure 7: Installation with sanitary adapters

4.7 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

5 Product specifications

5.1. Sensor specifications

	InPro 7250 ST/HT	InPro 7250 PFA
Measuring principle	inductive	inductive
Sensor data		
Measuring range	0...2000 mS/cm	0...2000 mS/cm
Measurement value deviation	ST: ± (0.5 % of meas. value + 25 µS) HT: ± (0.5% of meas. value + 1 µS)	± (0.5% of meas. value + 25 µS)
Cell factor nominal ¹⁾	2.175	2.30
Transfer ratio	120	120
Temperature sensor	Pt 1000	Pt 1000
Temperature response (t ₉₀)	approx. 5 min	approx. 5 min

Ambient conditions ²⁾

Process temperature (general purpose)	ST: -20 °C...+100 °C (-4 °F...+212 °F) HT: -20 °C...+180 °C (-4...+356 °F)	-20 °C...+125 °C (-4...+257 °F)
Process temperature (In Ex-classified area installation)	ST: - HT: -20 °C...+130 °C (-4...+266 °F)	-20 °C...+125 °C (-4...+257 °F)
Process pressure ³⁾	ST: max. 8 bar (max. 116 psi) HT: max. 20 bar (max. 290 psi)	max. 16 bar (max. 232 psi)

 **Note:** See also permissible pressure and temperature diagrams in appendix.

Sensor materials

Sensor (medium wetted) ⁴⁾	PEEK glass-fiber reinforced (GF30)	PFA not glass-filled
Sensor mounting thread (G 3/4")	stainless steel	stainless steel
O-ring	Viton®	FEP
Flat gasket	Viton	PTFE

Cable

Type	ST: single coaxial cable HT: double coaxial cable	simple coaxial cable
Cable jacket	ST: PVC HT: silicone	PVC
Available cable lengths	3, 5 and 10 m (9.8, 16.4 and 32.8 ft)	3, 5 and 10 m (9.8, 16.4 and 32.8 ft)

Approvals and certificates

ATEX, UKEX	ST: - HT: •	- •
FM	ST: - HT: •	- •
Quality/End control	•	•

¹⁾ This value depends on the installation.

The exact value must be determined with a calibration.

²⁾ The combination of high pressure, high temperature and/or aggressive process medium reduces the sensor life span.

³⁾ Depending on the installation.

⁴⁾ This material is not FDA listed.

6 Ordering information

6.1. Sensors

Sensor	Order no.
InPro 7250 ST/Pt1000/3 m	52 002 736
InPro 7250 ST/Pt1000/5 m	52 002 737
InPro 7250 ST/Pt1000/10 m	52 002 738
InPro 7250 HT/Pt1000/3 m	52 002 739
InPro 7250 HT/Pt1000/5 m	52 002 740
InPro 7250 HT/Pt1000/10 m	52 002 741
InPro 7250 PFA/Pt1000/3 m	52 005 423
InPro 7250 PFA/Pt1000/5 m	52 005 424
InPro 7250 PFA/Pt1000/10 m	52 005 425

6.2 Spare parts

Designation	Order no.
O-ring Viton	52 750 171
Flat gasket Viton	52 403 432
Hexagon nut (stainless steel)	52 403 433
O-ring FEP	52 750 192
Flat gasket PTFE	52 403 833

6.3 Process Connections and Accessories

– Flanges	Order Number
Flange DN 50/PN16	52 403 565
Flange ANSI 2"	52 403 567
Flange ANSI 3"	52 403 569
Flange DN50/PN16, PVDF, only for PFA version	52 403 946
Flange ANSI 2", incl. Sealing Plate PTFE	52 403 947
– Bushings	
Bushing R 1½"	52 403 446
Bushing R 1½", PVDF	52 403 447
Bushing R 2"	52 403 448
Bushing R 2", PVDF	52 403 449
Bushing 1½" NPT	52 403 450
Bushing 1½" NPT, PVDF	52 403 451
Bushing 2" NPT	52 403 452
Bushing 2" NPT, PVDF	52 403 453
– Sanitary Adapters	
Dairy adapter DN50	52 403 583
Aseptic adapter DN50	52 403 584
– InDip 550 Ind – Sensor holder	
InDip 550 ind PVC	52 403 579
InDip 550 ind PVDF	52 403 580
– Accessories	
Flat gasket (Viton)	52 403 432
O-ring (Viton)	52 750 171
Locknut (stainless steel)	52 403 433

6.4 Transmitter M700(X), module Cond Ind 7700(X)

Transmitters	Designation	Order Number
Transmitter base, ss (no modules)	M700S	52 121 174
Transmitter base, ss, Ex, VPW*, 100...230 VAC	M700XS/VPW	52 121 175
Transmitter base, ss, Ex, 24 VAC/DC	M700XS/24V	52 121 176
Transmitter base, coated (no modules)	M700C	52 121 171
Transmitter base, coated, Ex, VPW*, 100...230 VAC	M700XC/VPW	52 121 172
Transmitter base, coated, Ex, 24 VAC/DC	M700XC/24V	52 121 173

* VPW = VariPoWer

Conductivity		
measurement modules	Designation	Order Number
Conductivity (inductive) measurement module	Cond Ind 7700	52 121 186
Conductivity (inductive) measurement module, Ex	Cond Ind 7700X	52 121 187

6.5 Transmitter M400, (4-wire Transmitter)

Transmitter	Order Number
M400, Type 1 Cond Ind	52 121 495

6.6 Transmitter M400, (2-wire Transmitter)

Transmitter	Order Number
M400 2XH Cond Ind	30 256 307

Please contact your local METTLER TOLEDO representative for more information.

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US-Normen

Titel	Nummer	Ausgabe- datum
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

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Titel	Nummer	Ausgabe- datum
Explosive Atmospheres – Part 0: Equipment – General Requirements	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	CAN/CSA-C22.2 No. 61010-1	R2017

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Titre	Numéro	Date de publication
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

Normes canadiennes

Titre	Numéro	Date de publication
Explosive Atmospheres – Part 0: Equipment – General Requirements	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	CAN/CSA-C22.2 No. 61010-1	R2017

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