METTLER TOLEDO

DG113-SC

Glass electrode

Glaselektrode

Electrode de verre

Electrodo de vidrio

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Combined pH electrode

with movable sleeve junction for titrations in non-aqueous solutions.

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Care instructions

Attention

Due to transportation and storage the movable sleeve ring of the sleeve junction may stick to the glass and thus stop the flow of electrolyte. This may cause a very instable or drifting signal and may increase the response time.

Procedure

1. Close the filling aperture with the rubber cap.
2. Warm the movable sleeve ring with warm water and loosen it by hand by turning
3. Turn the electrode upside down and allow it to cool (loose sleeve junction).
4. Close the movable sleeve by a slight turning motion.
5. Check the fill level of the electrolyte and refill if necessary.
6. The electrode is ready to use.

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
1. Start-Up maintenance

- 1.1 The internal buffer solution (3) must cover the whole inside of the glass membrane (1). Remove air bubbles by gently shaking in vertical direction.
- 1.2 To equalize pressure remove rubber cap (7) or perforate its nipple (8).
- 1.3 Refill the reference electrolyte (6) (e.g. METTLER TOLEDO Nr. 51340052) to approx. 1 cm below the aperture (7).
- 1.4 To activate the glass membrane, place the electrode in deion. water for half a day.

2. Adjustment and calibration

The parameters zero point and slope are measured by adjustment. A new electrode should have parameter values according to 2.1–2.2. Instructions for the adjustment method are given in the titrator manual.

- 2.1 The zero point of the electrode (pH value at $E = 0 \text{ mV}$) is approx. pH 7. The electrode potential, measured in a buffer solution of pH 7, must be between -50 mV and $+50 \text{ mV}$.
- 2.2 The electrode slope (potential change per pH unit) at the pH range 4 to 9 must be greater than -55 mV/pH (at 25°C ($32\text{--}131^\circ\text{F}$)).
- 2.3 You should test the response time of the electrode. After 30 s the indicated electrode potential should not change more than 2 mV within the following 30 s.
- 2.4 If calibrating a used electrode the zero point should be between pH 6 and 8, the slope not less than -52 mV/pH , and the potential change after 60 s should not exceed 3 mV/30 s.

3. Storage and lifetime

- 3.1 Store the electrode in the reference electrolyte – the sleeve junction (4) has to be immersed – and close the aperture (perforated rubber cap is sufficient).
- 3.2 Never let the glass membrane dry out, otherwise it has to be reactivated.
- 3.3 The lifetime of the electrode is 1 to 3 years, depending on its use and treatment.

4. Maintenance instructions

- 4.1 The measuring range extends from pH 0 to 12, the tolerable temperature range from 0 to 60/100 °C. (32–140/212 °F).
- 4.2 The sample solution should never enter into the reference electrode through the ground glass junction. Therefore the level of the reference electrolyte must always be higher than the sample solution level.
- 4.3 Just blot the electrode after rinsing, never rub it dry.
- 4.4 Do not let the cable connector plug corrode.
- 4.5 A dried out glass membrane can be reactivated by placing it for several hours in 0.1mol/L HCl.

5. Cleaning procedures

- 5.1 Contamination of the reference electrolyte:
Drain the reference electrolyte and proceed according to 1.3 and 1.4
- 5.2 Other contaminations: place it for a few minutes in an ultrasonic bath with deion. water, ethanol or a mixed acid solution (e.g. aqua regia).

6. Regeneration

If the electrode does not meet the qualifications according to 2.4 anymore, try a corrosion of the glass membrane: place it 1 to 2 minutes in diluted HF solution (2% HF/5% HCl, e.g. METTLER TOLEDO No. 51340073), rinse thoroughly with deion. water and treat it as in 5.1.

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