Contents:

4 Membrane bodies InPro 5000
25 ml CO₂ electrolyte
1 O-ring set

Inhalt:

4 Membrankörper InPro 5000
25 ml CO₂-Elektrolyt
1 O-Ring-Set

Contenu:

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InPro® 5000 (i)
52 002 464

Mettler-Toledo AG, Process Analytics, Im Hackacker 15, CH-8902 Urdorf, Switzerland
Phone +41 44 729 62 11, Fax +41 44 729 66 36, Subject to technical changes.
06/2015 © Mettler-Toledo AG, Switzerland. Printed in Switzerland. www.mt.com/pro

Instructions for electrolyte and membrane exchange of the InPro® 5000 (i) CO₂ sensors

**Electrolyte exchange:**

After each calibration of the interior body (without membrane body) new electrolyte should be used.

**Membrane body exchange:**

We recommend to use a new membrane body for each batch (sterilization). Prior to membrane body exchange the interior body (pH electrode) must be calibrated.

When replacing the membrane body and electrolyte, strictly observe the following instructions:

1. Unscrew the cap sleeve from the shaft and carefully pull it from the sensor.
2. Usually the membrane body remains in the cap sleeve and should be ejected by pushing it with the flat finger tip. Otherwise pull the membrane body from the interior body (pH electrode).

Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.

3. Clean the interior body with electrolyte.
4. Check all O-rings for mechanical defects and replace if necessary.
5. Half-fill the membrane body with CO₂ electrolyte and make sure that all bubbles are removed. Air bubbles can be removed by gently knocking on the membrane body.
6. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and must be removed with a paper tissue.
7. Carefully slip the cap sleeve over the fitted membrane body and screw it down.
8. Remove expelled electrolyte with a paper tissue.
9. Sensor is now ready for installation, sterilization and process calibration.

Membrane Kit for InPro® 5000 (i) CO₂ Sensors
Membran-Kit für InPro® 5000 (i) CO₂ Sensoren
Trousse de rechange des capteurs CO₂ InPro® 5000 (i)

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