Automated Water Content Determination of Transformer Oils with InMotion KF

Introduction

Transformer oil or insulating oil is an oil that is stable at high temperatures and has excellent electrical insulating properties.

Limits of 30 to 35 ppm of water in the transformer oil are generally referenced, as water affects the dielectric properties of the insulation and the aging rate of the insulating materials. Low water content in the transformer oil is therefore a very important parameter for safe operation, reliability and aging of a transformer. In the extreme case, transformers can fail because of excessive water in the insulation.

Karl Fischer titration coupled with an oven sample changer is a well-established and accurate technique for the determination of water in transformer oils.

In this application a coulometric C3OS Karl Fischer titrator connected to an InMotion KF Pro was used to first perform a temperature scan of the different transformer oil samples in order to identify their ideal gas-phase extraction temperature and to subsequently measure their water content.

As air contains oxygen, which could react with organic samples, an inert gas (nitrogen) was used for the analysis of the transformer oil samples.

Sample preparation and procedures

- Scan principle: A sample with an unknown water gas-phase extraction temperature is heated with a constant heating rate from a defined start to a defined end temperature. The released water is monitored as function of temperature. The qualitative interpretation of the curve allows to identify the optimal gasphase extraction temperature.
- Gas-phase extraction principle: Nitrogen from a gas cylinder was used as the purge gas. The two-stage pressure regulator should be used so that the final pressure is in the range of 0.5 – 1 bar. The nitrogen is then passed by a gas stop valve through a drying unit filled with silica gel and subsequently through a drying unit filled with molecular sieve. The nitrogen will transfer the water into the titration cell.
- 3. Long Needle: Use the long needle that is available for the 10 mL vial size for the sample analysis. This long needle will allow purging through the lubricant and, thus, liberating all water from it.
- 4. *Blank preparation*: Apart from the water in the sample, the sample vial also contains atmospheric humidity; This makes blank value determination necessary.

Place three empty 10 mL blank vials into positions 1 - 3 of the InMotion KF and close them with a screw cap.



- Sample preparation: Add ~ 6.0 g of the transformer oil sample into a 10 mL vial and close it with the screw cap. Place all filled vials into the corresponding positions of the InMotion KF rack.
- 6. Start of analysis: Place an empty vial into the drift position on the rack. Start M746. The method will first perform a pretitration to remove excess water from the titration cell and will subsequently go into the Standby modus. Perform manual drift determinations until the online drift value is $< 5 \mu q/min$. Then press "Start sample". Alternatively you can also select "Automatic" as analysis start and define the "Drift stability" accordingly. The analysis will automatically start as soon as the online drift is below the defined value. The method will usually start with the following sequence: Drift - Blank - Sample. In most cases the drift is determined once before the sample loop and this drift value is used for all subsequent calculations.

Chemistry

Water, M = 18 g/mol.

 $\begin{array}{l} \mathsf{ROH} + \mathsf{SO}_2 + 3 \ \mathsf{RN} + \mathsf{I}_2 + \mathsf{H}_2\mathsf{O} = (\mathsf{RNH}) \bullet \mathsf{SO}_4 \ \mathsf{R} + \\ 2 \ (\mathsf{RNH})\mathsf{I} \end{array}$

Solutions

- Chemicals: 100 mL HYDRANAL[®] Coulomat AG Oven.
- Standard: 1 % oven water standard
- Sample: Transformer oils

Instruments and accessories

- Compact line titrator C30S, titration cell without diaphragm (30252662)
- Titration Excellence T7 (30252675) or T9 (30252676) with coulometric Karl-Fischer kit (30267113)
- XPE205 Analytical balance (30087653)
- InMotion KF Pro Oven Autosampler with 10 mL rack (30407502)
- Gas stop valve SV2 (30407442)
- LabX software
- Spatula

Results

Scan:

For best results, the gas-phase extraction temperature needs to be high enough to ensure fast and complete water release. However, if the measurement temperature is too high organic molecules are liberated by decomposition and react with the titrant (i.e. with iodine), ultimately leading to an overestimation of the water content. A major increase in drift, i.e. cumulative water content at high temperature is an indication of decomposition. For the investigated transformer oil samples decomposition started between 150 to 210 °C. Optimal gas-phase extraction temperature (oven temperature) was chosen approx. 30 °C below the decomposition temperature.

Water content determination:

Water content determination of the NYTRO transformer oil sample at an oven temperature of 170 °C and a measurement time of 600 s:

	R2. Water content	
	Blank: [ua] NYTPO Sample: [nnm]	
	Bidrik. [µg], NTIKO Sumple. [ppm]	
Blank 1	38.6	
Blank 2	35.7	
Blank 3	41.4	
1	8.67	
2	9.34	
3	8.43	
4	7.45	
5	8.62	
Mean	8.50 ppm	
S	0.68 ppm	
srel	8.013 %	

Oven temperature, measurement time and water content of three transformer oil samples:

	Oven T	Time	W	ater conte	nt
	[°C]	[s]	Mean [ppm]	s [ppm]	srel [%]
NYTRO	170	600	8.50	0.68	8.013
T22	120	600	23.27	0.83	3.561
T400	145	600	61.24	0.73	1.184

Remarks

Coulometric reagents have a limited water capacity. Water capacity of the HYDRANAL[®] Coulomat AG Oven is 700 mg H_2O per 100 mL reagent.

Waste disposal and safety measures

Dispose as organic solvent

References

http://www.mt.com/global/en/home/products/Labo ratory Analytics Browse/Product Family Browse titrators main/automated-titration-systems/KFoven-autosampler.html



Measured values



Scan NYTRO sample



Water content determination NYTRO sample 1/5

Consumption	Meas. value	H ₂ O	Drift	Time
mC	mV	μg	µg/min	S
0.0	340.3	0.0	0	0
64.7	334.2	6.0	182.2	1
300.3	290.9	28.0	846.1	2
522.6	6 177.2	48.8	1411.2	3
597.2	118.1	55.8	1297.6	4
616.7	118.7	57.6	679.7	5
637.7	124.6	59.5	420.5	6
663.2	118.4	61.9	334.3	7
694.6	3 112.5	64.8	304.6	8
717.0	99.7	66.9	283.2	9
717.1	101.2	66.9	155.2	10
718.8	104.9	67.1	60.9	11
732.3	105.3	68.4	25.4	12
1255.8	96.5	117.2	0	594
1255.8	97.5	117.2	0	595
1255.8	103.1	117.2	0	596
1257.3	101.7	117.4	0	597
1257.4	101.7	117.4	0	598
1257.4	101.6	117.4	0	599
1257.6	97.1	117.4	0	600
1257.6	99.9	117.4	0	600

NYTRO sample 1/5

Method

0

Titration Application Note

001 Title	
Туре	InMotion KF Coul.
Compatible with	C30S / T7 / T9
ID	M746
Title	NYTRO
Author	Administrator
Date/Time	11/07/2017 04:07:23 pm
Modified on	11/07/2017 04:07:23 pm
Modified by	Administrator
PIDIECI	No
50F	None
002 Drift determination	
Wait time	60 s
Duration	3 min
003 Sample (Blank)	
Sample	
Sample type	Blank
Number of IDs	1
ID 1	Blank
Entry type	Fixed pieces
Pieces	l pcs.
Weight per piece	1.0 g
	1.0
	25.0 C
	Wandar
004 Titration stand (InMotio	on KF)
Туре	InMotion KF
Titration stand	InMotion KF/1
Temperature ramp	No
Oven temperature	170 °C
Source for drift	Determination
	IO µg/min
Elow rato	
Transfer tube heating	No
005 Mix time	
Duration	60 s
000 Tituation (KE Ooul) [1]	
Type	Polarized
Sensor	DM143-SC
Unit	mV
Indication	Voltametric
Ipol	5.0 µA
Stir	
Speed	45 %
Control	
End point	100.0 mV
Control band	250.0 mV
Rate	Normal
Generator current	Automatic
I ermination	Dolay tires
Type Dolay time	Delay lime
Min time	10 S 600 s
Max. time	600 s
007 Calculation P1	
	InMotion blank value
Result unit	
Formula	R]=(ICEQ/10 712-
	TIME*DRIFT)/C

Constant C =

METTLER TOLEDO

Decimal places	1
Result limits	No
Record statistics	Yes
Extra statistical functions	No

008 End of sample

009 Record

Summary	Yes	
Results	No	
Raw results	No	
Resource data	No	
Method	No	
Series data	No	
0 Blank		

010 Blank

Name	Blank IMKF Coul
Value B=	Mean[R1]
Unit	μg
Limits	No

011 Sample

Sample	
Sample type	Sample
Number of IDs	1
ID 1	NYTRO
Entry type	Weight
Lower limit	0.0 g
Upper limit	5.0 g
Density	1.0 g/mL
Correction factor	1.0
Temperature	25.0 °C
Entry	Arbitrary

012 Titration stand (InMotion KF)

Туре	InMotion KF
Titration stand	InMotion KF/1
Temperature ramp	No
Oven temperature	170 °C
Source for drift	Determination
Max. start drift	10 µg/min
Carrier gas source	InMotion KF
Flow rate	80 mL/min
Transfer tube heating	No
013 Mix time	
Duration	60 s
014 Titration (KF Coul) [2]	
Sensor	
Type	Polarized
Sensor	DM143-SC
Unit	mV
Indication	Voltametric
logi	5.0 µA
Stir	
Speed	45 %
Control	
End point	100.0 mV
Control band	250.0 mV
Rate	Normal
Generator current	Automatic
Termination	
Туре	Delay time
Delay time	10 s
Min. time	600 s
Max. time	600 s

1

015 Calculation R2

Result	Cont. blank comp.
	(B in µg)
Result unit	ppm
Formula	R2=(ICEQ[2]/10.712-
	TIME[2]*DRIFT-B[Blank
	IMKF Coul])/(C*m)
Constant	C = 1
Decimal places	3
Result limits	No
Record statistics	Yes
Extra statistical functions	No

016 End of sample

017 Record

Summary	Yes
Results	No
Raw results	No
Resource data	No
Method	No
Series data	No

Method (Scan)

001 Title

Туре	Scan KF Coul.
Compatible with	C30S / T7 / T9
	Coop
	Scoll
Title	Scan NYTRO
Author	Administrator
Date/Time	04/25/2017 11:05:28 am
Modified on	04/25/2017 11:05:28 am
Modified by	Administrator
Protect	No
SOP	None
001	
002 Sample	
Sample	
Sample type	Sample
Number of IDs	1
	NYTRO
Entry type	Weight
Lower limit	0.0 g
Upper limit	5.0 g
Density	1.0 g/mL
Correction factor	1.0
Analysis start	Manual
Entry	Arbitrary
L () y	
003 Titration stand (InMotion KF))
Τνρε	InMotion KF
Titration stand	InMotion KE/1
Max start drift	
	ro µg/mm
Carrier gas source	
Flow rate	80 mL/min
Transfer tube heating	No
004 Mix time	
Duration	20 s
005 Titration (KE Coul) [1]	
Sensor	
Туре	Polarized
Temperature sensor	DM143-SC
Unit	mV
Indication	Voltametric
Inol	5.0
ipoi	υ.υ μη
Temperature program	
Start temperature	40 °C
Hoating rate	-0 0 2 °C/min
End temporature	2 0/11111
Ena lemperature	240 6
Stir	
Speed	35 %
opood	UU 10
Control	
End point	100.0 mV
Control band	250 mV
Generator current	automatic
006 Calculation R1	
Result	Sample size
Result unit	g
Formula	R1=m
Constant C=	1
Μ	M[None]
Z	z[None]
Desired places	4

No

Yes

No

007 Record

Summary	Per sample
Results	Per sample
Raw results	Per sample
Table of meas. values	Yes
Sample data	Per sample
Resource data	Per sample
H ₂ O - T	Yes
Drift - T	Yes
H ₂ O - T & Drift - T	Yes
Method	Yes
Series data	Yes
008 End of sample	
009 Record	
Summary	Yes

ourning	
Results	Yes
Raw results	Yes
Resource data	Yes
Method	Yes
Series data	Yes