METTLER TOLEDO

DX202-SC

Plastic reference electrode Kunststoff-Referenzelektrode Electrode de reference en matière synthétique

Plastic reference electrode

Plastic reference electrode suitable as an electrical reference in ion selective measurements as well as for pH or metal half-cell electrodes in aqueous media. Especially suitable for applications in hydrofluoric acid containing samples.

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
1. Changing the bridge electrolyte

Unscrew the outer shaft and pour the bridge electrolyte out. If a different electrolyte is added, ensure that there is no precipitation. It may be necessary to use intermediate solutions for the electrolyte change. Information regarding the selection of the bridge electrolyte can be found in the electrode specification bulletins or application brochures.

2. Storage and lifetime

Unscrew the plastic shaft and drain the bridge electrolyte. For storage screw on the plastic shaft and place electrode in the electrode sleeve filled with a small amount of deionized water. The service life of the reference electrode is shortened when the bridge electrolyte chamber constantly remains filled above the internal fiber junction. Ensure that the fiber junction is always kept damp.

3. Maintenance instructions

Keep outer fiber junction always clean and damp. If necessary clean the unscrewed plastic shaft with acid or base and rinse it thoroughly with deionized water. The junction can be cleaned with deionized water by applying a slight vacuum or small pressure via a rubber balloon.

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© Mettler-Toledo AG 2006 ME-51710469 Printed in Switzerland 0611/2.31
Mettler-Toledo AG