
METTLER TOLEDO

DS800-TwoPhase

Surfactant sensitive electrode for two-phase titrations in organic solvents

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METTLER TOLEDO

General information

The DS800-TwoPhase surfactant sensitive electrode is suitable for the indication of anionic and cationic surfactants either in pure aqueous solutions as well as for the titration according to the DIN standards EN14480, EN14668 und EN14669 in two-phase systems comprising organic solvents. The electrode is maintenance-free and resistant against solvents, which are recommended in the standards.

Specifications

Type of ion selective membrane:

Polymer membrane

Resistance of the membrane:

< 1 MΩ (in use), < 20 MΩ (dry)

Measurement range:

< 10⁻⁵ mol/L to critical micelle concentration (cmc)

Measurement temperature range:

0...50 °C

Optimum pH range:

1...12

Interferences:

In surfactant mixtures the most sorbing surfactant is the dominant species. Strong lipophilic anions and cations such as tetraphenylborate ions, ionic dyes can interfere.

Precision

For pure tensides ±2 mV, drift < 0.5 mV/min

Response times

When changing from lower to higher surfactant concentrations, the response time is approximately one minute; on change from higher to lower concentrations several minutes.

Materials

Electrode body of POM copolymer, membrane module of POM and solvent resistant sensor material, seals of silicon.

Electrolyte for DX200 reference electrode

3 mol/L KCl (for anionic and cationic surfactants).

Disposal

Dispose the metallic, plastic and glass parts separately.

Maintenance and service

Deposits on the surface of the membrane can be removed by thoroughly rinsing with deionized water. Obstinate deposits may be treated by short rinsing (5 – 10 sec) with pure ethanol without harming the module. After rinsing, conditioning in a dilute sodium chloride solution (0.01 mol/L) during approx. 10 minutes is required. If the sample medium is changed it may be advantageous to condition the electrode during 30 minutes in the new sample medium in order to achieve a stable potential. For short breaks in between measurements, the electrode should be kept in deionized water. For longer breaks, e.g. over night, it is recommended to store the electrode in a dilute sodium chloride solution (0.01 mol/L).

Membrane modules with a damaged membrane, excessive response times or unstable potentials (high membrane resistance) require the exchange of the complete electrode.

Storage

The DS800-TwoPhase electrode can be stored in air, preferably in the membrane protection cap included in the standard delivery. It is best stored in the supplied tube during longer periods. Before use the electrode is conditioned in dilute sodium chloride solution (0.01 mol/L) for at least 10 minutes. **Eventually shake the electrode once like a clinical thermometer.**

Reference electrode DX200: allow electrolyte to run out. For storage, close filling port with the rubber cap.

Solvent Behavior of DS800 Sensor Membrane		
solvent	swelling (weight gain)	
MIBK	48	%
CHCl ₃	600	%
Hexane	2	%
Acetone	93	%
Toluene	48	%
THF	dissolution %	

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© 03/2011 Mettler-Toledo AG, 51710531A
Printed in Switzerland

Mettler-Toledo AG, Analytical
Sonnenbergstrasse 74, CH-8603 Schwerzenbach, Switzerland
Tel. ++41 44 806 7711, Fax ++41 44 806 7350
Internet: www.mt.com

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