
English	Reference Manual QS3 Sample Changer for XPR Analytical Balances
Deutsch	Referenzhandbuch QS3 Probenwechsler für XPR-Analysenwaagen
Español	Manual de referencia Cambiador de muestra QS3 para balanzas analíticas XPR
Français	Manuel de référence Passeur d'échantillons QS3 pour balances d'analyse XPR
简体中文	参考手册 QS3自动进样器 适用于XPR分析天平
日本語	リファレンスマニュアル QS3サンプルチェンジャ XPR分析天 びん用

METTLER TOLEDO

Reference Manual **QS3 Sample Changer**

English

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Table of Contents

1	Introduction	3
1.1	Further documents and information	3
1.2	Explanation of conventions and symbols used	3
1.3	Acronyms and Abbreviations	4
1.4	Compliance information	4
2	Safety Information	4
2.1	Definitions of signal words and warning symbols	5
2.2	Product-specific safety notes	5
3	Design and Function	7
3.1	Function description	7
3.2	Overview	8
3.3	Components description	8
4	Installation and Putting into Operation	11
4.1	Selecting the location	11
4.2	Scope of delivery	12
4.3	Assembling the sample changer	12
4.4	Adjusting the sample changer	12
5	Operation	12
5.1	Accessing the weighing chamber	13
5.2	Loading the racks with sample vessels	13
5.3	Using other weighing methods	14
6	Maintenance	14
6.1	Cleaning	14
6.1.1	Cleaning the weighing pan and drip tray	15
6.1.2	Cleaning the tunnel, racks, and rings	15
7	Troubleshooting	16
7.1	Error symptoms	16
8	Technical Data	16
8.1	General data	16
8.2	Dimensions	17
	XPR analytical balance with full-height draft shield	17
	XPR analytical balance with half-height draft shield	18
9	Disposal	18
10	Accessories and Spare Parts	18
10.1	Accessories	18
10.2	Spare parts	20
	Index	21

1 Introduction

Thank you for choosing a METTLER TOLEDO instrument. The instrument combines high performance with ease of use.

1.1 Further documents and information

► www.mt.com/XPR-automatic

This document is available in other languages online.

► www.mt.com/QS3-RM

Instructions for cleaning a balance: "8 Steps to a Clean Balance"

► www.mt.com/lab-cleaning-guide

Search for software downloads

► www.mt.com/labweighing-software-download

Search for documents

► www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

1.2 Explanation of conventions and symbols used

Conventions and symbols

Key and/or button designations and display texts are shown in graphic or bold text, e.g., , **Edit**.

Note

For useful information about the product.

Refers to an external document.

Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.

1 Step 1

➡ Intermediate result

2 Step 2

➡ Result

1.3 Acronyms and Abbreviations

Original term	Explanation
EMC	Electromagnetic Compatibility
FCC	Federal Communications Commission
LPS	Limited Power Source
POM	Polyoxymethylene
RFID	Radio-frequency identification
RM	Reference Manual
sd	Standard deviation
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
UM	User Manual
USB	Universal Serial Bus

1.4 Compliance information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► <http://www.mt.com/ComplianceSearch>

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

► www.mt.com/contact

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Notice

2.2 Product-specific safety notes

Intended use

This dosing system is designed to be used in analytical laboratories by trained staff. The dosing system is intended for weighing and dosing powder or liquid samples.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Protective equipment

Chemical-resistant gloves

Goggles

Lab coat

Safety notes

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

WARNING

Injury and/or damage due to hazardous substances

Chemical, biological or radioactive hazards can be associated with the substances processed by the instrument. During dosing procedures, small amounts of the dosed substance may become airborne and penetrate the instrument or contaminate its surroundings.

The substance characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with the substance and take adequate safety measures, e.g., those stated on the safety data sheet provided by the manufacturer.
- 2 Make sure that every instrument part in contact with the substance will not get altered or damaged by the substance.

WARNING

Injury and/or damage due to reacting, flammable, or explosive substances

During the dosing procedure, substances could be combined and cause an exothermic reaction or explosion. This includes powders, liquids, and gases.

The sample characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with reacting, flammable, or explosive substances.
- 2 Ensure a working temperature low enough to prevent the formation of flames or an explosion.

CAUTION

Injury due to moving parts

- Do not reach into the working area while parts of the instrument are moving.

CAUTION

Injury due to sharp objects or broken glass

Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

NOTICE

Damage to the instrument

The instrument contains no user-serviceable parts.

- 1 Do not open the instrument.
- 2 In the event of problems, please contact a METTLER TOLEDO representative.

NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

3 Design and Function


3.1 Function description

The **QS3 sample changer** can be used in conjunction with the **Q3 dosing module** and can be installed on any **XPR analytical balance**. It contains 30 dosing positions (15 when used with sample vessels of diameter larger than 31.5 mm), and therefore provides increased throughput of your dosing activities. The sample changer can be used to dose liquids as well as powders. The dosing positions are covered by a tunnel draft shield with two handles, on the right and left side of the weighing chamber. The top door of the balance is replaced by a top panel especially designed to reduce the influence of air drafts on the weighing results.

Various adapter are available to accommodate a high variety of sample vessels. The vessels with adapters are placed on the racks, which are in turn installed on the rotating drive ring. Each dosing position is numbered for easy traceability. When dosing starts, the drive ring rotates and the dosing positions are brought one-by-one over the weighing pan for dosing.

Simple weighings, adjustments, and tests can still be performed when the sample changer is installed. To this purpose, a weighing pan cover and a weighing basket are provided with the instrument.

3.2 Overview


1	Side tunnel with handle, left	5	Side tunnel with handle, right
2	Home position	6	Top panel
3	Weighing pan sample changer	7	Back tunnel
4	Rack with 5 positions, 6 pcs		

3.3 Components description

Weighing pan and rail


A rail is fixed to the base ring. All adapters and sample vessels slide on this rail as the drive ring rotates. The rail has one gap above the weighing pan. When a sample vessel reaches this position, its weight hangs on the weighing pan, allowing weighing/dosing.

Drive ring


The drive ring is the rotating component of the sample changer. It contains six elongated slots to fit six racks, each with five dosing positions. The drive ring also presents a single round hole, the **Home position**, where simple weighing, adjustments, and tests can be performed. When the **Home position** is not used for weighing, cover the hole using the **Home position cover**.

The rotation of the drive ring is provided by a motor and a driving wheel, located at the rear of the sample changer.


Racks for sample vessels

The racks are used to hold the adapters (each containing a sample vessel). Up to six racks can fit on the sample changer, each with five dosing positions. Each position is numbered to allow dosing traceability.

Home position cover


The **Home position cover** is designed to avoid unwanted powder on the weighing pan. Whenever possible, cover the hole of the **Home position** with the **Home position cover**.

Weighing pan cover

When performing a simple weighing, adjustment, or test, you can use the weighing pan cover in the **Home position**.


If the drive ring rotates while the weighing pan cover or the weighing basket is installed, the sample changer or the load cell could be damaged. Always remove the weighing pan cover or the weighing basket as soon as the operation is completed. Place back the **Home position** cover.

Weighing basket


When using small weights that need to be positioned precisely in the middle of the weighing pan, use the weighing basket in the **Home position**.

If the drive ring rotates while the weighing pan cover or the weighing basket is installed, the sample changer or the load cell could be damaged. Always remove the weighing pan cover or the weighing basket as soon as the operation is completed. Place back the **Home position** cover.


Tunnel

To minimize the influence of air draft on the weighing, a tunnel is enclosing the weighing chamber and the drive ring holding the racks and sample vessels. The racks can be accessed by simply sliding tunnel on the right or left side of the weighing chamber.

Top panel

The top panel of the QS3 sample changer is designed to minimize the air draft in the weighing chamber.

4 Installation and Putting into Operation

This instrument must be installed by a METTLER TOLEDO service technician.

The QS3 sample changer is compatible with all XPR analytical balances with tall draft shield.

4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

Requirements of the location

Place indoors on stable table


Ensure sufficient spacing

Level the instrument

Provide adequate lighting


Avoid direct sunlight


Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Sufficient spacing for balances: > 15 cm all around the instrument

Take into account the environmental conditions. See "Technical Data".

4.2 Scope of delivery

QS3 sample changer

- Base ring with drive unit and slider rails
- Position pins, 2 pcs
- Base plate
- Sealing insert for side doors, 2 pcs
- Sealing insert for partition panel
- USB cable A – B
- Drive ring with numbered positions
- Rack with numbered positions, 6 pcs
- Top panel for sample changer
- Tunnel, 3 pcs
- Drip tray
- Weighing pan
- Weighing pan cover
- Weighing basket
- **Home position** cover
- Reference Manual

4.3 Assembling the sample changer

This instrument must be installed by a METTLER TOLEDO service technician.

For operation and maintenance, some parts of the sample changer can easily be disassembled, such as the tunnel, racks, and weighing pan.

See also

-  Loading the racks with sample vessels ► Page 13
-  Cleaning ► Page 14

4.4 Adjusting the sample changer

Adjusting the alignment – front to back

The alignment of the sample changer (front to back) is adjusted by the METTLER TOLEDO technician during installation.

Adjusting the alignment – left to right

The lateral alignment of the sample changer is adjusted by the METTLER TOLEDO technician during installation. If necessary, the alignment can be corrected, in particular when using asymmetric sample vessels, using the function **Adjust sample changer** on the balance terminal.

For further information, consult the Reference Manual (RM) of your XPR balance.

► www.mt.com/XPR-analytical-RM

5 Operation


For further information, consult the Reference Manual (RM) of your XPR balance.

► www.mt.com/XPR-analytical-RM

5.1 Accessing the weighing chamber

To access the weighing chamber, use the handle of the left or right tunnel to slide it towards the left or right, respectively.

5.2 Loading the racks with sample vessels

Each numbered position on the racks can be used for dosing. Use the appropriate adapter for your sample vessel to ensure stability. Do not use sample vessels without adapters. Do not leave empty adapters on the racks.

- 1 Use the handle to open the right or left side tunnel (1).
- 2 Remove the rack (2) accessible through the tunnel opening.
- 3 Unload the sample vessels from the rack, if applicable.
- 4 Load the rack with up to five sample vessels and adapters (if applicable).

 Note

For sample vessels of diameter larger than 31.5 mm, fill only every second position of the rack.

- 5 Install the rack at the right position on the drive ring.

 Note

To properly identify the sample vessels on the terminal, make sure that the numbering on the racks corresponds to the numbering on the drive ring.


 Note

For sample vessels of diameter larger than 31.5 mm, lift the adapter from the neighboring rack while loading the racks on the drive ring, if necessary.

- 6 Make sure that the rack is leveled.
- 7 Use the user interface to move the drive ring toward the left or right.
- 8 Repeat steps 2 to 7 until all the desired racks are installed in the sample changer.

 Note

If a sample position contains no sample vessel, the dosing sequence will skip this position and move to the next sample vessel.

For further information, consult the Reference Manual (RM) of your XPR balance.

► www.mt.com/XPR-analytical-RM

See also

 Accessories ► Page 18

5.3 Using other weighing methods

NOTICE

Damage to the instrument

The pin under the weighing pan cover and the weighing basket is inserted in the weighing pan during use. If the drive ring moves while the weighing pan cover or the weighing basket is installed, the sample changer or the balance can get damaged.

- 1 Remove the weighing pan cover or the weighing basket as soon as the weighing, test, or adjustment is complete.
- 2 Always check that the **Home position** cover is in place before moving the sample changer's position.

The sample changer is most often used with automated weighing/dosing. However, other weighing methods, test, and adjustments can still be performed with the balance.

- 1 Move the dosing module to its uppermost position.
- 2 Move the sample changer to the **Home position**.
- 3 Remove the **Home position** cover.
- 4  **WARNING: Damage to the instrument. Do not leave the weighing pan cover or the weighing basket on the sample changer while not in use.**
Install the weighing pan cover or the weighing basket.
- 5 You are free to use manual methods, tests, and adjustments.
- 6 Remove the weighing pan cover or the weighing basket.
- 7 Install the **Home position** cover.

For further information, consult the Reference Manual (RM) of your XPR balance.

► www.mt.com/XPR-analytical-RM

6 Maintenance

The appropriate maintenance interval depends on your standard operating procedure (SOP).

Please contact your METTLER TOLEDO representative for details about the available service options. Regular servicing by an authorized service technician ensures constant accuracy for years to come and prolongs the service life of your instrument.

6.1 Cleaning

NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.


For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► www.mt.com/lab-cleaning-guide

6.1.1 Cleaning the weighing pan and drip tray

- The racks, adapters, sample vessels, and **Home position** cover are removed from the drive ring.
 - One side tunnel is open and the front draft shield is removed.
- 1 Hold the drive ring (1) and lift it up above the weighing pan hooks.
 - 2 Remove the weighing pan (2).
 - 3 Remove the drip tray (3).
 - 4 Put down the drive ring.
 - 5 Clean the weighing pan and the drip tray.
 - 6 Reassemble all parts in reverse order.


6.1.2 Cleaning the tunnel, racks, and rings

For thorough cleaning, further components can be removed from the instrument. Clean the parts with a mild household cleaning agent, as required. Do not use cleaning agents containing abrasive components.


Disassembling the sample changer and cleaning

- The side tunnels are open and the front draft shield is removed.
- 1 Remove each side tunnel. Slightly spread the sides of the arch apart and lift.
 - 2 Remove the back tunnel.
 - 3 Remove the racks, the **Home position** cover, and the drive ring.
 - 4 If necessary, remove the weighing pan and the drip tray.
 - 5 Clean the tunnel, the drive ring, and the racks.
 - 6 If necessary, clean the weighing pan and the drip tray.
 - 7 Wipe the base ring.

Reassembling the sample changer

- 1 Install the drip tray.
- 2 Install the weighing pan.
- 3 Install the drive ring and the **Home position** cover.
- 4 If necessary, install the racks.
- 5 Install the back tunnel.
- 6 Install each side tunnel. Make sure that each tunnel wheel (1) is under the base ring.
- 7 Install the front draft shield.

7 Troubleshooting

For further information, consult the Reference Manual (RM) of your XPR balance.

► www.mt.com/XPR-analytical-RM

7.1 Error symptoms

Error symptom	Possible cause	Diagnostic	Remedy
The balance does not recognize the sample changer.	The USB connection of the sample changer was temporarily disconnected while the balance was running.	—	Verify that the sample changer is connected to the balance. Restart the balance by disconnecting the power cable and reconnecting after a few seconds.
The sample changer moves slowly and does not stop in the Home position .	—	—	Contact a METTLER TOLEDO representative.

8 Technical Data

8.1 General data

Weight (without packaging): 7030 g
Power consumption: 5 / 12 V DC $\pm$ 10%, 0.5 A

Protection and standards

Overvoltage category: II
Degree of pollution: 2
Range of application: Use only indoors in dry locations

Environmental conditions


Height above mean sea level: Up to 5000 m
Ambient temperature: +5 – +40 °C
Relative air humidity: 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

Storage conditions (in packaging)

Ambient temperature: -25 – +70 °C
Relative air humidity: 10 – 90%, non-condensing

8.2 Dimensions

Dimensions in mm.

XPR analytical balance with full-height draft shield

Dimensions	Vessel diameter ≤ 31.5 mm	Vessel diameter ≥ 33.5 mm
Minimum diameter of sample vessel opening	6 mm	6 mm
Maximum diameter of sample vessel	31.5 mm	52 mm
Minimum height of sample vessel ¹⁾	10 mm	10 mm
Maximum height of sample vessel ²⁾	140 mm	140 mm
Maximum number of sample vessels	30	15
Maximum load	6 kg	6 kg
Weight of QS3 sample changer with Q3 dosing module and XPR analytical balance with tall draft shield	16.4 kg	16.4 kg

¹⁾ when used with an adapter, height **excludes** the adapter

²⁾ when used with an adapter, height **includes** the adapter

XPR analytical balance with half-height draft shield

Dimensions	Vessel diameter ≤ 31.5 mm	Vessel diameter ≥ 33.5 mm
Minimum diameter of sample vessel opening	6 mm	6 mm
Maximum diameter of sample vessel	31.5 mm	52 mm
Minimum height of sample vessel ¹⁾	10 mm	10 mm
Maximum height of sample vessel ²⁾	100 mm	100 mm
Maximum number of sample vessels	30	15
Maximum load	6 kg	6 kg
Weight of QS3 sample changer with Q3 dosing module and XPR analytical balance with short draft shield	16.4 kg	16.4 kg

¹⁾ when used with an adapter, height **excludes** the adapter

²⁾ when used with an adapter, height **includes** the adapter

9 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

10 Accessories and Spare Parts

10.1 Accessories

Accessories are additional components that could help you in your workflow.

Please contact your METTLER TOLEDO sales representative for more details about using accessories with your product.

	Description	Order no.
Adapters		
	Racks for sample vessels with pins Set of 6 racks with positions numbered from 1 to 30 ø 24 mm	30246875

Vial adapters (POM), diameter < 15 cm

8.5 mm × 15 mm (5 pcs.)	30428901
9 mm × 25 mm (5 pcs.)	30428902
9.5 mm × 25 mm (5 pcs.)	30428903
10 mm × 25 mm (5 pcs.)	30428904
10.5 mm × 25 mm (5 pcs.)	30428905
11 mm × 30 mm (5 pcs.)	30428906
11.5 mm × 30 mm (5 pcs.)	30428907
12 mm × 20 mm (5 pcs.)	30428908
12.5 mm × 30 mm (5 pcs.)	30428909
13.5 mm × 20 mm (5 pcs.)	30428910
14.5 mm × 25 mm (5 pcs.)	30428911

Vial adapters (POM), diameter > 15 cm

15.5 mm × 25 mm (5 pcs.)	30428912
16.5 mm × 25 mm (5 pcs.)	30428913
17.5 mm × 25 mm (5 pcs.)	30428914
18.5 mm × 25 mm (5 pcs.)	30428915
19.5 mm × 25 mm (5 pcs.)	30428916
20.5 mm × 25 mm (5 pcs.)	30428917
21.5 mm × 25 mm (5 pcs.)	30428918
22.5 mm × 25 mm (5 pcs.)	30428919
23.5 mm × 25 mm (5 pcs.)	30428920
26 mm × 25 mm (5 pcs.)	30428926
27.5 mm × 25 mm (5 pcs.)	30428921
28.5 mm × 25 mm (5 pcs.)	30428922
29.5 mm × 35 mm (5 pcs.)	30428923
31.5 mm × 30 mm (5 pcs.)	30428924

Bottle adapters (POM)

33.5 mm × 15 mm (5 pcs.)	30459921
35.5 mm × 15 mm (5 pcs.)	30459922
37.5 mm × 15 mm (5 pcs.)	30459923
40.5 mm × 15 mm (5 pcs.)	30459924
44.5 mm × 15 mm (5 pcs.)	30459925
48 mm × 15 mm (5 pcs.)	30459926
52 mm × 15 mm (5 pcs.)	30459927

Tube adapters (POM)

Eppendorf 1.5 ml (5 pcs.)	30306209
Falcon 16.4 mm × 60 mm (5 pcs.)	30459992
Falcon 29.3 mm × 60 mm (5 pcs.)	30459991

Capsule adapters (stainless steel)

size 000 (5 pcs.)

30006416

size 00 (5 pcs.)

30006417

size 0 (5 pcs.)

30006418

size 1 (5 pcs.)

30006419

size 2 (5 pcs.)

30006430


size 3 (5 pcs.)

30006431

size 4 (5 pcs.)

30006432

10.2 Spare parts

	Order no.	Designation	Remarks
1	30542815	Tunnels, set	Including: 3 tunnels
2	11141772	Sample changer racks, set	Including: 6 racks with positions numbered 1–30
3	30460287	Leveling feet, set	Including: 2 leveling feet
4	30542823	Drip tray	—
5	30542816	Weighing pan sample changer	—
6	30542833	Weighing basket	—
7	30542817	Cover home position	—
8	30542818	Weighing pan cover large	—
9	30542819	Sealing insert top panel	Including: 5 sealing inserts
10	30542822	Panel top	—

Index

A

adjustment	14
altitude	16

C

compliance information	4
convention	3
cover	
home position	9
weighing pan	10, 14

D

dimension	17
disposal	18
draft shield	11

E

environmental condition	11, 16
-------------------------	--------

H

home position	
cover	9
humidity	16

I

install	
site	11

L

location	11
----------	----

M

manual weighing	14
-----------------	----

R

rail	8
------	---

S

safety information	5
sample changer	
adjust	12
simple weighing	14
storage	
condition	16
symbol	3
warning	5

T

technical data	16
----------------	----

temperature	16
test	14
top panel	11

W

warning symbol	5
weighing	14
weighing basket	10, 14
weighing pan	8
cover	10, 14
weighing basket	10

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© Mettler-Toledo GmbH 03/2021