METTLER **TOLEDO**

InPro 3200 (SG) pH Electrode

Instruction Manual

| | | |
|----|----------------------------|----|
| cs | Návod k obsluze | 2 |
| da | Brugsanvisning | 6 |
| de | Betriebsanleitung | 10 |
| en | Instruction manual | 14 |
| es | Manual de instrucciones | 18 |
| fi | Käyttöohje | 22 |
| fr | Instructions d'utilisation | 26 |
| hu | Használati utasítás | 30 |
| it | Istruzioni per l'uso | 34 |
| jp | 取扱説明書 | 38 |
| ko | 사용 매뉴얼 | 42 |
| nl | Gebruikershandleiding | 46 |
| pl | Instrukcja obsługi | 50 |
| pt | Manual de instruções | 54 |
| ru | Инструкция по эксплуатации | 58 |
| sv | Bruksanvisning | 62 |
| th | คู่มือการแนะนำ | 66 |
| zh | 说明书 | 70 |

InPro 3200 (SG)

pH elektroda

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Obsah

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
Instruction Manual

Contents

| | | |
|-----|---|----|
| 1 | Introduction | 15 |
| 2 | Safety instructions | 15 |
| 3 | Product description..... | 15 |
| 4 | Installation and commissioning | 16 |
| 5 | Operation..... | 16 |
| 5.1 | Calibration of the sensor and pH transmitter..... | 16 |
| 5.2 | Sterilization of sensor..... | 16 |
| 5.3 | Temperature sensor | 17 |
| 6 | Maintenance | 17 |
| 7 | Trouble-shooting | 17 |
| 8 | Disposal | 17 |
| 9 | Guarantee..... | 17 |

InPro is a registered trademark of Mettler-Toledo GmbH in Switzerland, the USA, the European Union and a further five countries.

1 Introduction

METTLER TOLEDO pH sensors type InPro® 3200 are low-maintenance electrodes with pressurized gel-electrolyte, silver-ion trap and built-in temperature sensor.

The InPro 3200 sensors are suitable for the same applications as the corresponding DPA(S) sensors without RTD, and provide an additional temperature signal for automatic temperature compensation of the pH signal during calibration and operation.

InPro 3200 sensors are autoclavable and in-situ sterilizable.

2 Safety instructions

The sensor consists partly of glass and should be handled carefully. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. Local safety regulations are applicable to the handling of the sensor. Please follow the guidelines enclosed for an installation in the Ex zone. Please take extra care during installation/fitting of the sensor to ensure that neither the polymer body/shaft of the sensor nor the thread (Pg 13.5) suffer any injury or damage.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

3 Product description

The inscription on each electrode contains the following information:

| | |
|-------------------------------|---|
| METTLER TOLEDO | Manufacturer of the electrode |
| InPro® 3200/SG/xxx/yyy | Type designation / SG = ancillary electrode / xxx = rod length in mm / yyy = temperature sensor |
| combination pH | Type of electrode |
| pH 0...14 | pH measurement range |
| 0...100 °C (140 °C) | Temperature range for operation (sterilization) |
| Order No. 52 00X XXX | Order number |
| Ex SEV14 ATEX 0168 X | Identification and number of certificate |
| Ex IECEx SEV 14.0025X | Identification and number of certificate |
| See instruction manual | Instruction (Observe the operating instructions) |

In addition, each electrode is allocated a serial number on the red part of the plug-in head to enable identification.

4 Installation and commissioning

1. On unpacking, check the sensor for possible mechanical damage of the pH sensitive glass, the diaphragm, the shaft and the VarioPin (VP) connector. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and briefly rinse the sensor with deionized water. Carefully remove the protection droplet of silicone on the diaphragm with the enclosed knife, as shown on the green pictures on the packing.
3. Again rinse the sensor well with deionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH sensitive glass, since this can lead to electrostatic charging and sluggish response times.
4. Carefully check the space behind the pH sensitive glass for the presence of any air bubbles and remove same by gently swinging the sensor in a vertical plane.
5. Please refer to the instruction manual of the relative housing in order to install the sensor correctly.
6. Connect the sensor to the pH transmitter, using a cable with VarioPin (VP) connector. Please refer to the wiring and terminal diagram accompanying the cable.
7. For electrodes with Pt 100 temperature sensors, 3-wire technique is recommended in order to obtain best measurement accuracy. Please refer to the cable and transmitter instruction manuals for the corresponding cable wire connections.
8. For InPro 3200 SG electrodes, with an additional Pt ring for solution ground, connect the corresponding blue wire in the VP-cable to your transmitter in accordance with instructions in the connection diagramm of VP-cable and/or in the transmitter instruction manual. For ORP measurement (sep. transmitter needed) connect solution ground with ORP input.

5 Operation

5.1 Calibration of the sensor and pH transmitter

For InPro 3200 sensors, a 2-point calibration is recommended. After removing the watering cap, the sensor is dipped in succession into two different buffer solutions with given pH reference values and the transmitter is calibrated to these buffer values.

5.2 Sterilization of sensor

In the case of sterile processes, such as in fermentation applications, the sensors are sterilized prior to the production cycle. Autoclaving of the whole fermenter/reactor unit including pH electrode is possible with this sensor type. Sterilization can also be effected in-situ by means of steam or superheated process medium in contact with the lower sensor parts.

5.3 Temperature sensor

The integrated Pt 100/1000 temperature sensor is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation, and **not** for the accurate indication or control of the process temperature.

6 Maintenance

1. Carefully rinse the sensor tip and diaphragm with deionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt™ electrolyte.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with deionized water or alcohol, and afterwards dry carefully.

7 Trouble-shooting

At sluggish response time (e.g. after a long and dry storage period), reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further sensor parts!) into Reconditioning (HF) Solution (article no. 51 319 001) for 5 min. Afterwards, rinse well with deionized water and let sensor tip and diaphragm soak in reference electrolyte for minimum 12 hours.

Rapidly increased resistance of reference electrode and/or drift could be signs of protein precipitation. Use Cleaning Solution Pepsin/HCl (article no. 51 340 068) to remove such precipitation. At repeated problems of this type, consider switching to sensors with Friscolyt™ based gel electrolyte.

8 Disposal

If no relative official regulations are in force, used or defective InPro 3200 sensors and the packaging can be disposed of as ordinary household waste.

9 Guarantee

On manufacturing defects, 12 months after delivery.

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Rev B. 08/2017 © Mettler-Toledo GmbH
Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, CH-8902 Urdorf