METTLER **TOLEDO**

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Instructions for Type 405-DPAS-SC pre-pressurized gel-reference pH Electrodes

1 Introduction

The METTLER TOLEDO Type 405-DPAS-SC combination pH electrode is a low-maintenance, pre-pressurized, autoclavable electrode with a gel-type reference. The reference element incorporates a silver-ion trap to prevent silver sulfide contamination of the reference junction. The gel-type reference is also available in a combination ORP and as a separate reference electrode.

2 Preparation

After unpacking, examine the electrode for any mechanical damage.

METTLER TOLEDO electrodes are factory tested and are shipped ready-to-use. (To maintain factory pressurization, the Type 405-DPAS-SC electrodes have a silicone seal covering the ceramic diaphragm; this must be removed with the enclosed blade prior to use.) Remove the hydration cap and remove any electrolyte crystals by rinsing the electrode with distilled water. Pat dry; do not rub as this may build up a static charge and lengthen the response time.

Eliminate any air bubbles around the internal element in the pH-sensitive membrane.

3 Operation

METTLER TOLEDO pH electrodes can be used with most pH/mV meters.

Follow manufacturer's instructions for use.

A two-point calibration is recommended, using fresh buffers of known pH values which bracket the sample pH; one buffer should be close to pH 7 to set the zero point. For fast routine analysis, a standardization using one buffer at a value near the sample pH can be done.

4 Sterilization

Type 405-DPAS-SC electrodes are steam-sterilizable in-situ and autoclavable. Follow operating instructions for the autoclave. It is advisable (but not necessary) to put the closure cap (p/n 00 201 1096) on the connector prior to autoclaving.

Repeated sterilization may cause the electrolyte to discolor slightly due to the elevated temperatures. This does not affect the electrode performance.

5 Maintenance

When properly cared for, METTLER TOLEDO electrodes should provide accurate results over thousands of uses.

Storage: Electrodes should be rinsed thoroughly after use and stored in 3M KCl or pH 4 buffer solution with added KCl. DO NOT store the electrode dry or in distilled or tap water; this dehydrates the pH membrane and the electrode must be rehydrated before use by immersing for several hours in 3M KCl.

Refilling: Type 405-DPAS-SC gel-reference electrodes are not refillable.

Protein Contamination: Clean the electrode with pH Electrode Cleaner for Proteins following instructions on the bottle.

Reactivating pH membrane: Use Reactivating Solution following directions CAREFULLY.

Servicing a clogged or dry reference junction: Immerse the electrode in a beaker of 3M KCl. Heat the solution to approximately 70 °C (160 °F) and allow electrode to soak in this solution as it cools to room temperature.

6 Troubleshooting

No response (symptom)

- Poor meter connection
 - Check connectors at the meter and at the electrode
- Cracked glass
 - Replace electrode

Slow response (symptom)

- Dehydrated pH membrane
 - Rehydrate (see Maintenance)
- Clogged reference junction
 - Soak in warm KCl (see Maintenance)
- Protein contamination
 - Use pH Electrode Cleaner (see Maintenance)
- Depleted gel layer
 - Use Reactivating Solution (see Maintenance)

Noise (symptom)

- Poor meter connection
 - Check connector
- Clogged reference junction
 - (See above)

Low slope (symptom)

- Dehydrated pH membrane
 - Rehydrate (see Maintenance)
- Contaminated membrane
 - (See above)

For further assistance, call our Customer Service Center at +1-800-352-8763.

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Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, CH-8902 Urdorf