Contents:

1. Product description
- 1.1 Equipment and scope of delivery
- 1.2 Drawing and parts list
2. Use
- 2.1 2-point dry calibration
- 2.2 Off-line calibration in open laboratory vessels
3. Conditions of warranty

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1.1 Scope of delivery: CaliCap with two plastic inserts, this instruction manual.

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METTLER TOLEDO Market Organizations:

BR Mettler-Toledo Ind. e Com. Ltda, Tel. (011) 4166 74 00
CH Mettler-Toledo (Schweiz) GmbH, Tel. 044 944 47 60
D Mettler-Toledo GmbH, Prozeßanalytik, Tel. (0641) 507 444
F Mettler-Toledo Analyse Indust. S.A.S, Tel. (01) 47 37 06 00
USA Mettler-Toledo Process Analytical Inc., Tel. (781) 301 8800

CaliCap – Calibration Tool for InPro 8000 Turbidity Sensors
52 800 248

For more addresses of METTLER TOLEDO Market Organizations please go to: www.mt.com/pro-MOs


12/2016 © Mettler-Toledo GmbH, Switzerland

Mettler-Toledo GmbH, Process Analytics
Im Hackacker 15, 8902 Urdorf, Switzerland
Tel. +41 44 729 62 11, Fax +41 44 729 66 36 www.mt.com/pro

ENGLISH

CaliCap – Calibration Tool for InPro 8000 Turbidity Sensors

1.2 Drawing of initial assembly with parts list

1 Stainless steel holder DIN1.4305
2 Screws (a, b)
3 Plastic insert black polypropylene
4 Plastic insert, light-colored polypropylene

InPro is a registered trademark of the Mettler-Toledo Group in Switzerland and a further twelve countries.

CaliCap is a trademark of the Mettler-Toledo Group.

2. Use

CaliCap™ is a precision calibration tool for METTLER TOLEDO scattered light sensors of the InPro® 8000 Series operating in combination with the METTLER TOLEDO turbidity transmitter M800 process 1-channel or Trb 8300. This tool is also suitable for use with the previous product series of FSC402 transmitters and OFS12 sensors. CaliCap serves to carry out a reproducible 2-point check of METTLER TOLEDO turbidity sensors without the use of liquid calibration solutions but is also suitable for calibrations with reference solutions in open laboratory vessels (beakers). The fault tolerance for a dry check comes within the limits of $\pm 5\%$ of the measuring value. Depending on plant SOP, it is recommended to carry out a reproducible 2-point dry check before each new production batch, or every 3 months for continuous production.

2.1 Dry check

- 1.) Loosen the screws **2a** and **2b** but do not remove them completely.
- 2.) Fully withdraw both plastic inserts **[3]** and **[4]** and put the light-colored insert **[4]** into where the black one **[3]** has been beforehand. Retighten screw **2b**.
- 3.) From the other side, push the sensor into the stainless steel holder **[1]** until it comes to rest against the light-colored plastic insert **[4]**. Fix the sensor with screw **2a**.

This constellation is intended to simulate a high turbidity value, whereby light is reflected by the surface of the light-colored insert.

- 4.) Vary the signal amplification (slope value) through the manual calibration function of the transmitter until a measuring value of 100 % appears on the display (see chapter "Manual Calibration" in the Instruction Manual for the transmitter.) Note the new value for the signal amplification.
- 5.) Release screw **2b** and replace the light-colored insert **[4]** by the black one **[3]**. Retighten screw **2b**. Leave the sensor in the same position.

This constellation is intended to simulate the zero point. Vary the offset value manually until a measuring value of "0" is shown on the display. Note the new offset value of the transmitter.

- 6.) Remove the sensor from the CaliCap and reassemble the tool as it was in its originally supplied condition.

The new values noted for signal amplification and for offset can later be used during dry checks when carrying out routine checks of the complete measuring system consisting of transmitter, sensor, fixed cable and possibly extension cable.

⚠ Any changes within the measuring system itself, such as replacement of the transmitter, sensor or cable, the use of a new optical fiber, an extension or even of a different CaliCap unit make it necessary to carry out a new calibration of the complete system.

2.2 Off-line calibration in open laboratory vessels

⚠ Prior to use of the CaliCap and InPro 8000 sensors in liquid samples, a check is to be made to ensure material compatibility of all wetted parts, or otherwise to abstain altogether from operation in aggressive media.

For off-line calibration, the stainless steel holder **[1]** together with only the black insert **[3]** is needed when dipping the sensor into the reference solutions. Due to the optimized design of the black plastic insert, any disturbing light reflections from the sides or bottom of the laboratory vessel are eliminated, thereby speeding up the calibration procedure.

The solutions/media should be stirred during calibration to prevent any sedimentation of solid particles.

Carry out steps 1, 2, 3 and 5 as described in **Section 2.1**. Steps 2 and 3 guarantee a reproducible constellation of sensor and CaliCap. Afterwards, the sensor with the CaliCap in place is to be dipped into all calibration solutions consecutively and calibration carried out as described in the Instruction Manual of the transmitter (see chapter "Manual Calibration" or "Multipoint Calibration"). Between the individual dippings into the different solutions, the sensor including the fitted CaliCap holder is to be rinsed each time with clear water or solvent.

When calibration has been completed, rinse the CaliCap and the sensor thoroughly with clear water and then dismount the sensor from the CaliCap unit completely. Rinse, clean and dry the individual parts **[1], 2a, 2b, 2c and [3]** before reassembling the CaliCap into its initial condition.

3. Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for the period of one year from date of delivery ex works.

If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, carriage paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiry of the period of warranty, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

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