METTLER TOLEDO

**Turbidity Sensors
Trübungssensoren
Sondes de turbidité**

InPro® 8600 i Series

InPro 8600 i Series
30 113 529

English

page 3

Deutsch

page 29

Français

page 55

Turbidity Sensors

InPro 8600 i Series

Instruction manual

© It is forbidden to reprint this instruction manual in whole or part.

No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other recordings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH-8902 Urdorf.

All rights reserved, in particular reproduction, translation and patenting/registration.

Content

1	Introduction	6
2	Safety	7
2.1	Introduction	7
2.2	EC type examination in accordance with directive 97/23/EC	8
2.3	Material Specification	9
2.4	EHEDG Certificate of compliance	10
2.5	Intended use	11
2.6	Inappropriate use	11
2.7	Basic principles	11
2.8	Responsibilities, organizational measures	11
2.8.1	Responsibilities of operator	11
2.8.2	Responsibilities of personnel	12
2.9	Product-specific hazards	12
2.9.1	Manipulation of and maintenance work on the sensor	12
2.9.2	Installation in pressurized systems	13
2.10	Residual hazards	13
2.10.1	Rupture of connections	13
2.10.2	Heat protection	13
2.10.3	External influences	13
2.11	Emergency measures	13
2.12	Safety measures	14
2.13	Modifications	14
3	Product description	15
3.1	Scope of delivery	15
3.2	Packing	15
3.3	Sensor designation	15
3.4	Product overview	16
3.5	Description of function	17
3.5.1	Measuring principle	17
3.5.2	Construction	17
4	Installation, operation and maintenance	18
4.1	Sensor installation	18
4.2	Electrical connections	20
4.3	Operation	21
4.3.1	Process conditions	21
4.4	Maintenance	22
4.4.1	Cleaning the sensor head	22
4.4.2	Calibration	22
4.4.3	Checking with an inspection unit	22
4.4.4	Tuchenhangen-Varivent® process adapter check and change of O-rings	22
4.4.5	Measurement window check	22
4.4.6	Humidity check	22
4.5	Storage	23
4.6	Disposal	23
5	Product specifications	24
6	Conditions of warranty	26
7	Ordering information	27

How to use this instruction manual

This instruction manual is an integral part of the METTLER TOLEDO InPro 8600i sensors and contains notes and instructions that are important for safety and operation.

All persons working on or with the InPro 8600i sensor must have first read and understood the sections appropriate to the work in hand.

Please read this instruction manual carefully before using the InPro 8600i sensor. Keep this document close to the unit, so that operating personnel may easily be able to refer to it at any time.

Please first read section 1 "Introduction" and section 2 "Safety".

Proprietary designations

The following are proprietary names and, for the sake of simplicity, will be mentioned in this instruction manual without the registration marking, e.g.®:

- InPro® is a registered trade mark of Mettler-Toledo GmbH, CH-8606 Greifensee, Switzerland.

Explanation of sensor designations

The general designation InPro 8600i used in this instruction manual refers to the following turbidity sensor models:

- **InPro 8600i/D1:** digital interface version, single-angle
- **InPro 8600i/D3:** digital interface version, dual-angle plus colour

1 Introduction

Thank you for buying **the InPro 8600i Series sensor from METTLER TOLEDO**. The construction of the InPro 8600i Series sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment. **Therefore, the instruction manual must be read and understood by the persons involved before work is started with the sensor.**

In addition to this instruction manual please also note the following:

- all local safety regulations
- all instructions and warning remarks in the publications of the products that are used in conjunction with the sensor
- all safety precautions for the plant into which the sensor will be installed.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Safety

2.1 Introduction

This instruction manual contains all the information needed for safe and proper use of the InPro 8600i Series sensor. The instruction manual is intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

The instruction manual must be stored where it is constantly accessible and available to any person working with the InPro 8600i.

Warning notices and symbols

Warning notices and informations are marked with the following symbols:

Danger! Warning of a danger that can lead to extensive material damage, to death or grave bodily injury.

Caution! Warning of a possible dangerous situation that can lead to light bodily harm and/or material damage.

Notice: Information referring to technical requirements. Non-adherence can lead to malfunctions, uneconomic working and possibly also to loss of productivity.

2.2 EC type examination in accordance with directive 97/23/EC

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2.4 EHEDG Certificate of compliance

2.5 Intended use

METTLER TOLEDO InPro 8600i Series sensors are intended solely for the measurement in liquids in industrial applications.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed
- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions
- observance of local legislation
- consultation with METTLER TOLEDO Process Analytics in the event of any uncertainties.

2.6 Inappropriate use

Any utilization other than the above mentioned, as well as any utilization with the technical data, is taken as being not in conformance with the regulations. The operator bears the sole risk for any damage caused by such utilization.

2.7 Basic principles

- The sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

Notice: Manipulation of the sensors may only take place after it has been ensured that no process medium can escape through the sensor in the event of incorrect manipulation. For this reason, the complete system must be emptied and vented in advance (**safe condition**).

2.8 Responsibilities, organizational measures

2.8.1 Responsibilities of operator

- The operator is under obligation only to permit persons to work with InPro 8600i sensors, who are familiar with the basic requirements of work safety and accident prevention, and who have been instructed in the handling of the sensor. This instruction manual serves as the basic document.
- In addition to the instruction manual there are also generally valid legal and other binding regulations for work safety and accident prevention as well as for environmental protection, and these must be provided by the operator and instructed to personnel using the sensors.

- The operator/user must be fully aware of safety and potential danger factors during work with the sensor and this awareness must be checked by the operator at regular intervals.
- Measures must be taken to ensure that the sensor is only operated in a safe and fully functional condition.

Caution! Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

2.8.2 Responsibilities of personnel

- All persons whose duty it is to operate the sensor are under obligation to read chapter 1 "Introduction" and chapter 2 "Safety" as well as the warning notices in this instruction manual.
- In addition to the instruction manual, generally valid legal and other binding regulations for work safety and accident prevention must be adhered to.
- Any method of working which is doubtful from a safety perspective and which exceeds the operation according to regulations must be omitted.

Notice: Before every start-up, the sensor must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- perfect functioning
- authorization for use in conjunction with other plant equipment and resources.

Caution! A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

2.9 Product-specific hazards

2.9.1 Manipulation of and maintenance work on the sensor

Notice: Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (depressurize, no explosion risk, empty, rinse, vent, etc.).

It is principally necessary to wear personal protective outfit such as protective goggles and clothing.

Only such maintenance work as is specifically described in this instruction manual may be performed on the sensor.

Only original spare parts from METTLER TOLEDO may be used for replacing defective components (see chapter 7).

Danger! Non-compliance with the prescribed maintenance instructions can endanger personnel and the environment.

2.9.2 Installation in pressurized systems

Notice: The maximum permissible temperature and pressure specifications must not be exceeded. The relative specifications are given on the type plate of the sensor.

2.10 Residual hazards

2.10.1 Rupture of connections

Sensor connections can become loosened through the effects of vibration.

Notice: All sensor interfaces including process adapter, and the interfaces between optical window and sensor head, need to be checked periodically by the customer and kept functional.

Danger! Leaky connections can allow process medium, cleaning solution or control (compressed) air to escape to the environment, presenting a hazard for persons and the environment.

2.10.2 Heat protection

Caution! The sensor is not equipped with heat protection. During steam-sterilization procedures in pipes, the surface of the sensor can reach high temperatures and cause burns.

2.10.3 External influences

Notice: Objects falling on the sensor can damage or destroy the unit, or cause leaks etc. Lateral forces may damage or destroy the unit.

2.11 Emergency measures

Notice: Always observe and comply with local regulations!

2.12 Safety measures

Notice: Always observe and comply with local laws and regulations!
These are not an integral part of this instruction manual.

The operator is responsible for the instruction of personnel. Additional copies of this instruction manual can be ordered from the equipment supplier. As an integral part of the sensor, this instruction manual must at all times be readily accessible to users at the point of operation of the sensor.

The operator is obliged to inform the supplier/manufacturer of the sensor immediately about any safety-relevant incidents, or observations made, during use of the sensor.

Danger! Incorrect manipulation and/or instruction errors can lead to potential hazards for persons and for the environment.

Before every start-up the sensor must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- defective cables and lines etc.
- authorization for use of the sensor in conjunction with the associated plant resources.

Caution! A defective sensor must never be installed or put into operation. Poor containment, leaky connections etc. or non-compliant installation of the sensor can lead to escape of process medium and thereby to a potential threat to life.

2.13 Modifications

Notice: No attachments or modifications to the sensor housings are allowed.

Danger! The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments and alterations or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

3 Product description

3.1 Scope of delivery

Standard supply of the sensor is made up of the following:

- instruction manual
- certificates in accordance with the specifications
- blanking plate

On receipt of the shipment, check immediately:

- The sensor and accessories for any sign of transport damage. Report any damage immediately to the carrier and to your supplier.
- The type designation on the sensor body.
- For completeness of the supply. Please notify your supplier immediately if the shipment is incomplete or in any way incorrect.

Notice: A defective sensor must neither be installed nor put into service.

3.2 Packing

The packing consists of cardboard and foam plastic. Keep the packing for possible later use, such as for storage or transport of the sensor. However, if you wish to dispose of the packing, please observe your local regulations in this respect.

3.3 Sensor designation

Type designation

A sensor can be identified by the article number and the model name on the label located on the sensor.

InPro 8600 i/xy

x: D = digital connection of the sensor with a transmitter M800 1-channel

y: 1 = single-angle version (25°)


3 = dual-angle plus colour version (25° and 90° plus colour)

Serial number

The serial number can be found on the label on the sensor.

Sensor certifications

The label located on the sensor provides the information about the sensor certifications.


Easy sensor identification due to the label

Legend:


- 1 Designation
- 2 Article number
- 3 Material specifications
- 4 Max. permissible pressure
- 5 Max. permissible temperature
- 6 Serial number
- 7 Sensor certification
- 8 Information about disposal

3.4 Product overview

The InPro 8600i/D* sensors can only be used in combination with the METTLER TOLEDO transmitter M800 1-channel.

Dimension drawing InPro 8600i/D*

3.5 Description of function

3.5.1 Measuring principle


The InPro 8600i series of sensors are used for **the determination of undissolved particles or substances in solutions**. They are designed for measuring low to medium levels of turbidity and colouration.

The sensors in the InPro 8600i series operate on **the scattered-light principle** in order to determine turbidity. Undissolved particles or substances located in the beam path of the sensor scatter the impinging light. Detectors arranged at an angle of 25° (single-angle version) or 25° and 90° (dual-angle version) measure the scattered light. The ratio of scattered light to transmitted light is output as a turbidity value. The turbidity value calculated in this manner is proportional to the concentration of undissolved particles or substances.

The InPro 8600i/D3 also has a colour measurement function (with or without turbidity compensation). Absorption of blue light is measured for the colour measurement. The graduation of yellow is determined according to the Mitteleuropäische Brautechnische Analysenkommission (MEBAK) directives. This allows to measure turbidity and graduation of yellow with one instrument.

3.5.2 Construction

The sensor consists of a sensor head and a sensor body. The sensor head, which is in contact with the medium, is attached to the pipe using a Tuchenhagen Varivent® process adapter. The instrument electronics are contained in the sensor body.

Assembly of the InPro 8600i sensor

4 Installation, operation and maintenance


Caution! In all installation work described below, ensure that the equipment in which the sensor is to be installed is in a non-hazardous condition (depressurized, empty, rinsed, vented, etc.).


4.1 Sensor installation

The sensor is equipped with a Tuchenhagen-Varivent® process adapter, which is used to attach the sensor to the process pipe in dimensions from DN40 to DN 150. Prior to the installation, the sealing surfaces (sensor and pipe adapter) need to be cleaned and checked for damages. Further, it needs to be verified that no obstacles are present within the installation cavity of the sensor.

The installation of the sensor should be in a **vertical standpipe**. If the sensor will be installed in a horizontal pipe, it must be in a horizontal position and the process pipe has to be filled completely with liquid during the measurement.


Installation on a vertical line (preferred)


Installation on a horizontal line

There is a mark on the sensor head which shows the **flow direction of the media**. The sensor needs to be installed in this direction.

Mark for the installation into a pipe

The sensor must be fitted to the pipe **at least 2 meters (6.6 ft) away from stray light sources** (such as sight glasses).

Notice: All of the illustrated items must be fitted when the sensor is installed to ensure proper operation and splash protection (IP65). Also ensure that cable feedthroughs are firmly secured and the cables cannot move. Seal unused cable feedthroughs with blanking plugs.


4.2 Electrical connections

Caution! Some components inside the sensor are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero-potential before any operation on the wiring terminal. Be careful during preparing of operation, handling and operation of the sensor. This work is to be accomplished in principle only by qualified technical personnel.

Before opening the cover, it must be ensured that the equipment, in which the sensor is installed, is in a safe condition (depressurized, empty, rinse, vent, etc.).

Pin assignments

InPro 8600i/D1/D3

InPro 8600i/D*

Terminal assignment	Terminal	Assignment	Cable lead
KL2	7	shield	—
	8	RS485 B (-)	green
	9	RS485 A (+)	yellow
KL1	10	24V DC –	white
	11	24V DC +	brown
	12	shield	green/yellow

Notice: The InPro8600i/D* turbidity sensor operates with a 24 VDC power supply only. When only an AC power supply is available on site, the InPro8600i/D* turbidity sensor may connect to the 24 VDC terminals from the M800 1-channel transmitter. However, when 24 VDC power supply is available on site, the InPro 8600i/D* should be directly connected to the 24 VDC power source.

M800 1-channel

Terminal assignment	Terminal	Function	Cable lead
TB3	2	GND5V	green/yellow
	3	RS485B	green
	4	RS485A	yellow
	7	24V	brown
	8	GND24V	white

In order to ensure optimum measurement results the following points must be considered:

- Run the connection cable through the protection.
- Connect the wires to the sensor.
- Hand-tighten the cable protection.

4.3 Operation

In case of any failure or malfunction, the equipment in which the sensor is installed, needs to be in a safe condition (depressurized, empty, rinse, vent, etc.).

4.3.1 Process conditions

- The plant pressure should never exceed the specification of the delivered sensor.
- The process temperature should never exceed the specification of the delivered sensor.

Caution! Exceeding the specified maximum pressure and/or the specified maximum temperature represent a very high safety risk. Please read the additional safety instructions before installation.

- Avoid air and gas bubbles inside the sensor, they cause disturbances. Noise and drift of the measurement would be the result (air bubbles are not expected at pressures higher than 2 bar in aqueous solutions).
- The sensor must be air cooled if the process temperature is persistently higher than 70° C (158 °F). The air connection ports must be used for this purpose (available as optional accessories).

4.4 Maintenance

4.4.1 Cleaning the sensor head

The measurement principle provides a large degree of compensation for contamination of the sensor head. However, the contamination can become so severe, depending on the operating conditions and media, that the compensation is no longer adequate and the sensor head must be cleaned.

1. Place the system or equipment in a safe state (depressurised, empty, rinsed, vented or ventilated, etc.)
2. Disconnect electrical power
3. Dismount the sensor
4. Clean the sensor head using a soft, lint-free cloth and a mild cleaning agent
5. Install the sensor

4.4.2 Calibration

The METTLER TOLEDO InPro 8600 i/D* turbidity sensor is factory calibrated prior to shipment. The factory calibration is done with 12-points formazin calibration over the entire measuring range. The factory calibration data is stored in the sensor and does not have to be entered manually. The LED light source in the InPro 8600 i is very stable and there is no change in intensity over time. The ratio measuring principle compensates for sample color change, or window fouling. Hence, it is not necessary to re-calibrate the InPro 8600 i/D* turbidity sensor. When re-calibration with formazin is required, the InPro 8600 i/D* turbidity sensor should be returned to the factory. Please consult your local METTLER TOLEDO supplier for assistance.

In-situ one point process calibration (slope and offset) are possible for the InPro 8600 i/D* turbidity sensor. Please refer to the METTLER TOLEDO M800 1-channel transmitter instruction manual, section "Turbidity Calibration" for details.

4.4.3 Checking with an inspection unit

The sensor can be checked using an inspection unit with a built-in solid reference (turbiglass) every month or when necessary. The procedure is described in the corresponding instruction manual.

4.4.4 Tuchenhagen-Varivent® process adapter check and change of O-rings

The O-ring of the process adapter needs to be checked monthly and during every service call. METTLER TOLEDO recommends that medium-wetted seals should under all circumstances be replaced every 6 months.

4.4.5 Measurement window check

The sapphire windows need to be checked monthly and during every service call. Take special attention to damage, contamination and sealing.

4.4.6 Humidity check

The sensor actively monitors the humidity level in its electronics compartment. The humidity level can be read out at the M800 transmitter. If the humidity level in the sensor increases, the sensor will issue an alarm. Under these circumstances the sensor's readings might be unreliable.

Removing moisture from the sensor's electronics compartment necessitates servicing. Please contact METTLER TOLEDO for assistance.

4.5 Storage

Please inspect the instrument immediately after receiving for eventual transport damages. The sensor should be stored in a dry, clean and protected area until time of installation.

Storage for a reinstallation

If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

Shipment of the sensors

Please clean the instrument carefully before shipment (e.g. for revision/repair). Please use a fixed packaging to protect the instrument against transport damages. In ideal case the original packaging should be used.

4.6 Disposal

The disposal is to the responsibility of the user. Sensors contain electronic components that require special disposal without any health hazard for humans, and with no risk of harm for the environment.

Notice: Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

5 Product specifications

Specifications InPro 8600i

Turbidity

Measuring principle	Scattered light, red light InPro 8600i/D1: 25° (single-angle version) InPro 8600i/D3: 25°/90° (dual-angle version)
Wavelength	650 nm (light-emitting diode [LED])
Measuring range	0 ... 400 FTU
Accuracy	0.01 FTU (measuring value < 1 FTU) 1 % of the measured value (meas. value > 1 FTU)
Repeatability	0.01 FTU
Resolution	0.01 FTU
Response time (T90)	< 2 sec
Units	FTU, NTU, EBC, ASBC, ppm SiO ₂ , mg/L SiO ₂
Factory calibration	12-point calibration over the entire measuring range, based on Formazin standards

Colour (InPro 8600i/D3)

Measuring principle	Absorption, blue light
Wavelength	460 nm (light-emitting diode [LED])
Measuring range	0 – 30 EBC
Measurement accuracy	± 2 % of the measured value
Repeatability	0.01 EBC
Resolution	0.01 EBC
Response time (T90)	3 sec at EBC change > 10
Units	EBC
Factory calibration	6-point calibration

Process conditions

Admissible pressure range	max. 16 bar (232 psi)
Temperature (process media)	InPro 8600i/D1: 0 ... 100 °C (32 ... 212 °F) Max. 120 °C (248 °F), short time, 1 h
	InPro 8600i/D3: 0 ... 80 °C (32 ... 176 °F), at 25 °C (77 °F) ambient temperature with air cooling 0 ... 70 °C (32 ... 158 °F), at 25 °C (77 °F) ambient temperature without air cooling Max. 120 °C (248 °F), short time, 1 h
Admissible temperature range (ambient)	0 ... 60 °C (32 ... 140 °F)

Materials and dimensions

Sensor head (wetted)	stainless steel (1.4404/316L)
Wetted measurement window	sapphire
Sensor housing	stainless steel (1.4404/316L))
O-ring (flange wetted)	EPDM
Sensor connection	in-line adaption (VARINLINE® or compatible Tuchenhagen N 50/40)
Surface finish	Ra ≤ 0.8 µm
Protection rating	IP65
Cable lengths	5, 10 and 20 m (16.4, 32.8 and 65.6 ft) (RS485 communication)
Pipe diameter (process)	DN40 ... DN150

InPro 8600i/D*

Communication with M800 1-channel	digital (RS485)
Power through M800 1-channel	24 V DC (± 2 V), reverse-polarity protected
Sensor power consumption	max. 300 mA

Certificates

Quality control and final inspection	•
PED	•
EHEDG	•
CE	•

6 Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for the period of one year from date of delivery ex works. If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, carriage paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiry of the period of warranty, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

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METTLER TOLEDO Market Organizations

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Great Britain

Mettler-Toledo LTD
64 Boston Road, Beaumont Leys
GB -Leicester LE4 1AW
Phone +44 116 235 7070
e-mail enquire.mtuk@mt.com

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Norway

Mettler-Toledo AS
Ulvenveien 92B
NO-0581 Oslo Norway
Phone +47 22 30 44 90
e-mail info.mtn@mt.com

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Vietnam

Mettler-Toledo (Vietnam) LLC
29A Hoang Hoa Tham Street, Ward 6
Binh Thanh District
Ho Chi Minh City, Vietnam
Phone +84 8 35515924
e-mail MT-VN.CustomerSupport@mt.com

Management System
certified according to
ISO 9001 / ISO 14001

Subject to technical changes.
© Mettler-Toledo GmbH, Process Analytics
07/16 Printed in Switzerland. 30 113 529

Mettler-Toledo GmbH, Process Analytics
Industrie Nord, CH-8902 Urdorf, Switzerland
Phone +41 44 729 62 11, Fax +41 44 729 66 36

www.mt.com/pro