
Pre-Operation Treatment of Electrodes

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


Pre-operation treatment of electrodes

- a When unpacking, check electrode for possible mechanical damage to the membrane (1), the diaphragm or the body. Electrodes subject to complaints must be returned to us or to our agent in the original box indicating the invoice number.
- b Take off watering cap (7) and flush the electrode tip with distilled water. After rinsing with distilled water, wipe the electrode with a tissue. Do not rub as an increase of the response time might result.
- c Remove rubber cap which covers the electrolyte filling aperture (4) or cut off nipple (T) to ensure pressure equalization with the atmosphere.
- d Eliminate any air bubbles inside the membrane space (6) of the measuring electrode by smoothly swinging the unit in the vertical plane.
- e To renew a dried out glass membrane, soak it for several hours in the appropriate storage solution. It is recommended to store glass electrodes in dist. water, combined and reference electrodes in the reference electrolyte solution (as under h).
- f After a long and dry storage period, the glass membrane can be reactivated with hydrofluoric acid (order no. 51 340 073): The electrode membrane is dipped for 1 min. in this solution followed by rinsing with distilled water and soaking for 12 hours in the storage solution.
- g Electrodes with ground sleeve diaphragm: Slightly lift the ground sleeve so that some reference electrolyte drops out.
- h Check level of electrolyte (5) which should be approx. 10 mm below refilling aperture (4).
Refill, if necessary with:
 - 3M KCl + AgCl (order no. 51 340 045) for electrodes with Ag/AgCl wire
 - 3M KCl (order no. 51 340 049) or Viscolyt™ (order no. 51 340 235) for high temperature electrodes with Argenthal cartridge
- i Electrodes with gellified electrolyte without refilling port are treated as per §§ a, b, d, e, f above. For storage we recommend a solution 3 mol/l KCl.
- k Before calibrating or taking measurements, attention should be paid to the instructions of the relative pH instruments manufacturer. Additionally, please refer to the METTLER TOLEDO brochure "Practice and Theory of pH Measurement".

Viscolyt is a trademark of the METTLER TOLEDO Group.

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All rights reserved