
METTLER TOLEDO

InFit® 764-50
InFlow™ 764-20/-51

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Instruction manual

Contents

1.	Introduction	5
2.	Important notes	5
2.1	General.....	5
2.2	Safety precautions.....	5
3.	Description of product.....	6
4.	Installation and start up procedures	6
4.1	Installation	6
4.1.1	Mounting in reactors and pipes	6
4.1.2	Installing and removing the housing	7
4.1.3	Installing and removing the electrode	8
4.1.4	Checking for correct installation	9
4.2	Start up.....	10
4.2.1	Calibrating the electrode assembly.....	10
4.2.2	Pressure compensation	10
5.	Routine operation	11
6.	Maintenance	11
7.	Trouble shooting	12
8.	Specifications.....	13
8.1	Delivery	13
8.2	Technical specification.....	13
8.3	Spare parts and accessories	14
8.3.1	Housing.....	15
8.3.2	Electrodes	16
8.4	Warranty.....	16
9.	Appendix.....	17
A.	Drawing InFit® 764-50	17
B.	Drawing InFlow™ 764-20/-51	18
C.	Instruction for welding-in the threaded sockets	19

1. Introduction

These operating instructions describe how to use the insertion housing **InFit® 764-50**, and the flow-through housing **InFlow™ 764-20/-51**. In the case of non-standard models, additional information is to be found on the attached specification sheet.

The housings are state of the art, highgrade engineering products, tested by METTLER TOLEDO. Nevertheless, improper handling can be dangerous.

Attention! Work on the housing should be assigned only to trained personnel. Observe local regulations concerning the safety of people and property. Be sure to follow the instructions in this operating manual carefully.

2. Important notes

2.1 General

Immediately on receipt, check that the housing is complete and in good condition. Notify your supplier of any damage or deficiency.

Please also refer to your supplier for further information on ordering spare parts and accessories.

2.2 Safety precautions

- InFit® 764-50 and InFlow™ 764-20/-51 housings are designed only to contain METTLER TOLEDO screwcap electrodes. Any other kind of use could be dangerous and is not permitted.
- The materials used for the housing are described in Section 8 (and also on the specification sheet for nonstandard versions). Make sure the materials are suitable for the required application.
- To ensure that the housing is correctly installed and maintained, follow the instructions given in this manual. Incorrect handling of the housing can result in a broken electrode and leakage from the vessel or piping.
- Before doing anything to the installed housing, ensure that the process facility is in a safe condition (release pressure, empty, rinse, vent, purge, etc.).
- Use only clean electrodes, housings and weld-in sockets. Replace damaged seals and housing components.
- Before starting up, always check the measuring system. Inspect the housing/electrode assembly and check for leaks from housing and apparatus.
- If ever in doubt, consult your supplier.

3. Description of product

The insertion housing **InFit® 764-50** and flow-through housing **InFlow™ 764-20/-51** (the latter composed of an insertion housing InFit® 764-50 and a flow-through chamber 764-20 (DN25) serve as enclosures for METTLER TOLEDO screwcap electrodes used for industrial pH and redox measurements. The housings can be sterilized *in situ* and may be mounted obliquely or vertically in reactors or pipes. Electrodes with screw caps are easy to change and allow the electrode cables to be used repeatedly. The large electrolyte reservoir reduces maintenance to a minimum.

All parts of the housing in contact with the process are of stainless steel DIN 1.4435 (other material options are available on request); all remaining parts are of nickel plated brass (ring nut and cylinder are optionally available in stainless steel). The top piece is made of polypropylene. Silicone rubber O-rings (temperature resistant up to 150 °C) seal in the process medium; all other seals are of neoprene. The flow-through chamber is fabricated from stainless steel DIN 1.4408.

4. Installation and start up procedures

4.1 Installation

For the installation please observe the information in the following chapters.

4.1.1 Mounting in reactors and pipes

Choose a position for the housing that allows sufficient access. Take account of the housing's overall length (dimensions, see Appendix A and B).

The insertion housing InFit® 764-50 is attached to a weld-in socket (25 mm inside diameter). On reactors the sockets can be on the side or on top. The weld-in sockets can be used for pipes larger than 50 mm nom. diameter. For pipes of smaller diameter, the flow-through chamber 764-20 should be used. The chamber is attached to the two connecting pieces provided. These connecting pieces are welded to the pipe.

When mounting the sockets and the flow-through chamber, make sure the angle of the housing to the horizontal is **greater than 15°**; otherwise electrolyte will run out of the electrode reservoir. This will falsify measurements and corrode the housing. Also, be sure to use sockets that are long enough to allow the ring nut on the housing to be tightened.

The following lengths are recommended:

Weld-in sockets, straight	
Reactor wall thickness (mm)	socket length (mm)
6 – 13	40
13 – 16	50
16 – 25	60

Weld-in sockets 15° angle	
Reactor wall thickness (mm)	socket length (mm)
3 – 8	40
8 – 16	48
16 – 22	55

Brief instructions on welding are contained in Appendix C.

Attention! Observe relevant safety regulations for welding.

4.1.2 Installing and removing the housing

Installing:

1. Before installing, clean the immersion tube and the weld-in sockets.
2. Insert the housing in the sockets and tighten the ring nut by hand.
Important! The seal between housing and socket is an O-ring. **Therefore hand tighten only; do not use any tools.**

Removing:

1. **Before removing the housing, ensure that the process facility is in a safe condition (release pressure, empty, rinse, vent, purge, etc.).**
2. Undo the ring nut holding the housing to the weld-in socket. Then remove the housing.

4.1.3 Installing and removing the electrode

Attention! If plastic adapter "A" and metal ring nut "B" are loosened or tightened in the wrong order, the electrode may break.


Figure 1 "Installing and removing the electrode"

Removing:

- 1. Before removing the electrode, ensure that the process facility is in a safe condition (release pressure, empty, rinse, vent, purge, etc.).** The compensation pressure can be relieved either by slightly loosening the valve, or interrupting and venting the pressure supply.
- 2. First detach the connector plug from the electrode.** To do this, pull off the cable grommet "D", unscrew sleeve "C" and disconnect the cable connector from the electrode.
- 3. When removing the electrode or opening the housing, always undo the plastic adapter "A" first, then the metal ring nut "B".** Then remove the electrode.

Installing:

1. Remove the watering cap from the electrode tip and rinse with water.
2. Remove the elastic band and stopper from the electrolyte filler hole. The elastic band must taken off completely; the electrode must not be installed with the elastic band attached. Then insert the electrode into the housing.
Important! Note the **printed mark on the electrode “Position electrode this side up”** and turn the electrode to the correct position.
3. Always tighten the metal ring nut “B” first, then the plastic adapter “A” (ring nut “B” must be tightened only when plastic adapter “A” is loose).

Important! Hand tighten only, do not use any tools.

4. Finally, attach the plug connector.
To do this, push the jack of the cable from above through the sleeve “C” and screw the jack to the socket on the electrode. Then screw the sleeve fingeright to the plastic adapter “A”. Mount the split grommet “D” on the cable and press the grommet into the sleeve.

Important! The sleeve protects the connection between electrode and cable against dirt and moisture, and also against mechanical forces. **The protective sleeve should therefore always be used.**

Details on connecting the cable to the pH/redox transmitter are to be found in the operating instructions for the transmitter.

Important! Having fitted the electrode, never invert the housing.

4.1.4 **Checking for correct installation**

Important! Each time before starting up, check the measuring system. Inspect the electrode assembly and examine for leaks from housing and apparatus (also see Section 7).
Do not begin operation until the measuring system has been checked and any necessary corrective action taken.

4.2 Start up

4.2.1 Calibrating the electrode assembly

To calibrate the electrode assembly, first detach the housing (see “4.1.2 Installing and removing the housing”).

There is no need to remove the electrode from the housing when calibrating, nor does pressure compensation have to be turned off. The electrode and housing can be dipped directly in the buffer solutions.

Details of the calibration procedure are given in the operating instructions for the electrode and the transmitter.

Important! After installing the assembly, check it for leaks.

4.2.2 Pressure compensation

Important! The pressure in the reference electrode must always be **higher (by 0.2 to 2 bar)** than that of the medium (reactor contents or pressure in pipe) to ensure an outward flow of the electrolyte.

Remember that the hydrostatic pressure of the vessel must be taken into account and added to the reactor pressure. Pressure compensation may therefore be necessary even in a non-presurized reactor.

The pressure compensation is provided by the air pump or some other oil-free pressure source. If the latter is chosen, replace the valve assembly with the pressure connector assembly provided. The pressure can be read from the gauge.

Important! If the air pump is used for this purpose, the pressure must be frequently checked and adjusted.

The compensation pressure may be reduced by slightly loosening the valve assembly or by interrupting and venting the pressure supply.

5. Routine operation

The electrolyte in the reference electrode diminishes steadily as it flows out through the diaphragm. The level must therefore be checked regularly.

Top-up the electrolyte when it is close to the bottom of the reservoir bulb. To do this, remove the top piece of the housing and remove the electrode (see "4.1.3 Installing and removing the electrode"). For further details on routine use of the electrode, consult the electrode operating instructions.

Important! If pressure is provided by the air pump, the pressure must be frequently checked and adjusted.

Details on using the transmitter are contained in the operating instructions for the transmitter.

6. Maintenance

The electrode, housing and weld-in socket must be kept clean. Replace any damaged seals or components.

The insertion housing InFit® 764-50 can be sterilized in situ with the electrode installed.

Attention! Autoclaving or otherwise heating of the complete assembly is not allowed.

Details on maintenance of the electrode are contained in the electrode operating instructions.

7. Trouble shooting

A defective or poorly installed housing will leak. Testing for leaks is done with the electrode in place. For this, set the pressure in the housing to 6 bar, using the air pump. A leakproof housing loses less than 0.5 bar in 10 hours.

To check the seal between socket and housing, the reactor/pipe must be air-pressurized. Any air escaping can be detected by means of a leak-detecting spray.

Attention! Do not do anything to the housing before the pressure in the reactor/pipe and in the housing has been released again.

A leaky seal can be cured by cleaning the surfaces and lubricating the seals with a silicone-free grease. Faulty seals must be replaced.

Important! All seals are in the form of O-rings or resilient gaskets, so tightening hard with the aid of tools is not necessary.

Attention! Tightening the metal ring nut afterwards may cause the electrode to break (the ring nut must not be tightened unless plastic adapter "A" is loose (see "4.1.3 Installing the electrode").

Damaged or faulty parts of the housing must be replaced for safety reasons. In the case of non-standard versions, refer to the additional information in the appended specification sheet, and contact your supplier.

8. Specifications

8.1 Delivery

A **standard delivery** comprises the following items:

- Housing InFit® 764-50 or
InFlow™ 764-20/-51
- Protective sleeve with cable grommet 5 mm
- Bicycle pump
- Electrolyte top-up syringe
- Pressure connector assembly
- Cable grommet 7 mm
- Operating instructions
- Specification sheet for nonstandard versions

Options:

- Protective cage

8.2 Technical specification

Pressure range: 0 to 6 bar

Temp. range: Immersion tube: -30 to 130 °C,
steam-sterilizable

Top piece: -30 to 80 °C

Mounting: Weld-in socket internal dia. 25 mm

Thread: G1 1/4", 1 1/4" NPSM

Flow-through chamber DN25

Materials: (for non-standard versions, also see
specification sheet)

Top piece:

Nickel-plated brass, polypropylene,
glass, Teflon®

Wetted parts:

Stainless steel DIN 1.4435

Surface finish N6

silicone O-rings

Flow-through chamber DN 25:

Cast stainless steel DIN 1.4408

Dimensions: see Appendix A and B

Weights: InFit® 764-50: ~2 kg
Flow-through chamber DN 25: ~1 kg

8.3 Spare parts and accessories

To order spare parts and accessories, contact your supplier.


Figure 2 "Exploded drawing"

8.3.1 Housing

Important! With non-standard versions, see additional information on the specification sheet. Use only spare parts shown on the sheet. Use standard spare parts if the specification sheet does not indicate otherwise.

Note: Key to numbered parts shown in the exploded drawing.

Spare parts for standard version:

No.	Description	Material	Order No.
10	Teflon saddle	Teflon®	00 764 1053
20	Metal cylinder	nickel-plated brass	00 764 1054
30	Ring nut M60	nickel-plated brass	00 764 1058
40	Glass cylinder	glass	00 764 1064
50	Gasket 764-32.2	neoprene (CR)	00 764 1065
60 60/10	Top piece complete Plastic Adapter	polypropylene polypropylene	00 764 2074 00 764 2073
70	Gauge 0 - 6 bar compl.		00 764 2045
70/10	Valve body 764	nickel-plated brass	00 764 1062
70/20	Valve nut	nickel-plated brass	00 764 1061
70/30	Nut	nickel-plated brass	00 764 1060
70/40	Gauge 0 - 6 bar		20 100 1001
70/50	Gauge gasket	copper	20 102 1000
70/60	Valve		20 102 1001
70/70	O-ring 6.40 x 1.90 mm	nitrile (NBR)	20 300 1002
80	Lower probe H=70 w. fixed cage	steel 1.4435	00 764 2011
80	Lower probe H=70 no cage	steel 1.4435	00 764 2012
80	Lower probe H=100 no cage	steel 1.4435	00 764 2013
80	Lower probe H=150 no cage	steel 1.4435	00 764 2014
80	Lower probe H=200 no cage	steel 1.4435	00 764 2015
80	Lower probe H=250 no cage	steel 1.4435	00 764 2016
80	Lower probe H=300 no cage	steel 1.4435	00 764 2017
80	Lower probe H=350 no cage	steel 1.4435	00 764 2018
80	Lower probe H=400 no cage	steel 1.4435	00 764 2019
80	Lower probe H=450 no cage	steel 1.4435	00 764 2020
80	Lower probe H=500 no cage	steel 1.4435	00 764 2021
80/40	O-ring 10.77 x 2.62 mm	silicone (MVQ)	20 301 1132
80/50	O-ring 18.64 x 3.53 mm	silicone (MVQ)	20 301 1003
91	Cable grommet 5 mm	neoprene (CR)	00 764 1236
92	Sleeve	polypropylene	00 764 1283
110	Pressure connector		20 102 2000
170	Cable grommet 7 mm	neoprene (CR)	00 764 1237
180	Flow-through chamber DN 25	steel 1.4408	00 764 1087
190	Gasket	klingerite	00 764 1353
200	Connecting piece (welding)	steel 1.4435	00 764 1089
210	Nut St 18/8/2.5	steel 1.4435	00 764 1091

Optional Accessories:

No.	Description	Material	Order No.
20	Metal cylinder	steel 1.4435	00 764 1055
30	Ring nut M60	steel 1.4435	00 764 1059
A	Protection cage ¹⁾	steel 1.4435	00 764 1046
	Pt100 for flow-through chamber 5m		10 100 3104
	Pt100 for flow-through chamber EX 5m		10 100 3105

1) available for lower probes with lengths "H" = 100 mm or longer.

8.3.2 Electrodes

Suitable for use are pH electrodes of type 465-50-S7 and redox electrodes of type Pt4865-50-S7.

The length of the electrode for a given housing is defined as follows:

Electrode length = H + 50 mm ("H": see Appendix A and B)

For further details on electrodes, see documentation on the electrode or ask the supplier of your housing.

8.4 Warranty

The housings are of high technical quality and undergo a policy of continuous design review to incorporate the latest advances. Their reliability is ensured by a thorough final inspection prior to leaving our factory.

The warranty is valid for one year from the date of delivery and covers any defects due to faulty materials or manufacture.

Not covered by the warranty are normal wear and tear and any damage caused by improper use (e.g. chemical incompatibility of the materials, etc.).

The warranty extends only to replacement or repair of deficient products, at our discretion.


The warranty is void if the customer or others modify in any way the products supplied by us.

Defects must be reported to the supplier immediately upon discovery, and in all events within the warranty period.

9. Appendix


A. Drawing InFit® 764-50

Dimensions in mm

B. Drawing InFlow™ 764-20/-51

Dimensions in mm

C. Instruction for welding-in the threaded sockets

The threaded sockets are shipped ready for mounting, the bore -H7²⁾ being finish-ground.


It is due to the material that the inner diameter of the threaded sockets may slightly change when the adaptor is welded into the counter-piece (reactor wall or top plate). To prevent the bore of the threaded socket from becoming too small or too large we recommend make a small welding seam ($a = \text{max. } 3 \text{ mm}$) in one operation. You can obtain the best results with the TIG¹⁾ welding process. After the welding-in procedure the diameter of the bore needs to be checked using a gauge -H7²⁾ and, if necessary, slightly reamed out with a finishing hand-reamer.

1) TIG = Tungsten inert gas (TIG welding)


2) Tolerances -H7:

$\varnothing 18\text{-}30 \text{ mm: } +0.021 \text{ mm}$
 -0.000 mm





Weld-in socket

$\varnothing 25 \text{ mm straight}$

$\varnothing 25 \text{ mm angular}$

Mounting position, straight*into reactor top plate**into tube horizontal***Mounting position, angular***into reactor wall or tube**into reactor top plate or tube*

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Subject to technical changes. 01/15 © Mettler-Toledo AG. Printed in Switzerland. 59 906 035
- CH** **Mettler-Toledo (Schweiz) GmbH**
- CN** **Mettler-Toledo Instruments (Shanghai) Co. Ltd.**
- D** **Mettler-Toledo GmbH**
- F** **Mettler-Toledo, Analyse Industrielle S.A.S.**
- UK** **Mettler-Toledo LTD**