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Uso conforme

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Operating instructions

Intended Use

METTLER-TOLEDO InLab® pH electrodes are intended for precise pH measurement in the laboratory and in the field at atmospheric pressure.

Technical Data

The exact specifications for the pH and temperature range are printed on the electrode shaft and can be found at www.mt.com/electrode-guide.

General Safety Information

When working with chemicals, all relevant safety regulations provided by the manufacturer and the laboratory must be observed.

Product-specific Safety Information

The electrode is only intended for the purpose specified above. The manufacturer does not assume liability for any damage that results from use outside of the intended purpose.

The electrode is made partly of glass; there is therefore a risk of injury in the event of breakage. Failure to observe the handling instructions below may result in distorted results and electrode damage. Do not scratch the glass membrane, rub it dry with a cloth, or touch it with your fingers. Do not use the electrode to stir the measurement medium or touch the bottom or walls of the measuring vessel.

Further Information

Further information is available at www.mt.com/electrode-guide.

Manufacturer

Mettler-Toledo GmbH
Im Langacher 44, 8606 Greifensee, Switzerland
Phone +41 22 567 53 22, Fax +41 22 567 53 23
Internet www.mt.com

Commissioning

The electrode is supplied ready for use. When unpacking, please carry out a visual inspection to check that the glass is intact. Use a suitable cable to connect the electrode to the pH meter (see Fig. 1 and 2). Now remove the watering cap (see Fig. 3) and remove the silicone seal with the knife supplied (see Fig 4). Only remove the seal when you intend to use the electrode. Rinse the electrode with distilled water. If there is an air bubble inside the electrode tip, eliminate it by shaking the electrode up and down (as you would a clinical thermometer). The electrode is then ready for calibration and measurement.

Operating instructions

Storage

After use, rinse the electrode well with distilled water. Replace the watering cap (see Fig. 3), filled with electrolyte or InLab storage solution (order no. 30111142), (the glass membrane must be wetted). The electrode should be stored upright and at room temperature.

Troubleshooting

Make sure that ...

- all cables are connected correctly and the connections are clean.
- there is no air bubble inside the electrode tip.
- fresh buffer solutions are used for calibration.

For further information, visit

www.mt.com/electrode-guide -> pH TroubleShooter

Note: The electrode must be recalibrated after every treatment with a cleaning solution. After dry storage, the electrode must be wetted in electrolyte for approx. 24 hours.

Disposal

In accordance with the requirements of the European Directive 2002/96/EC on waste electrical and electronic appliances (WEEE), this product must not be disposed of as household refuse. This also applies equally for countries outside of the EU in accordance with their applicable national regulations. Please dispose of this product in accordance with the local regulations in a separate container for electrical and electronic devices. If you have questions regarding this, please contact your responsible authority or the dealer from which you purchased the product. If you pass the product on to a third party (e.g., for private or commercial/industrial use), then this regulation must also be passed on accordingly. Thank you for contributing to the care of the environment.

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