

# GEP Recommendation

**Company** ABZ pharma  
**Contact** Roger Silber  
**Department** Lab head  
**Site** Bern, Switzerland  
**City / Zip** bern

Date: 4/19/2023

## User Requirements

<b>Measuring parameter</b>	pH	<b>Expected pH value</b>	7
<b>Measuring location</b>	Changing location (Portable)	<b>Sample consistency</b>	Liquid
<b>Impact of inaccurate measuring</b>	Medium	<b>Sample composition</b>	Aqueous
<b>Measurement tolerance</b>	= ±0.05pH	<b>Expected sample temperature</b>	25 °C
<b>Measurement storage</b>	Via PC transfer	<b>Sample volume</b>	30 ml
<b>Number of samples per day</b>	1		

## System Recommendation

<b>Meter</b>	Seven2Go Pro S8-Meter
<b>Sensor</b>	InLab®Routine Go-ISM
<b>Accessory</b>	uGo carrying case
	License EasyDirect pH
<b>Calibration standards</b>	NIST/DIN Buffer pH 6.865
	NIST/DIN Buffer pH 9.180
<b>Verification standard</b>	Technical Buffer pH 7.00
<b>Maintenance solution</b>	Storage Solution InLab



## Services Recommendation

Services by supplier	
<b>Installation and Qualification</b>	Starter Pac
<b>User training</b>	Expert with certification
<b>Preventive maintenance</b>	Comprehensive Care
<b>Calibration &amp; adjustment</b>	Yearly

Routine operation by user - pH Sensor		
<b>Calibration standard</b> <sup>1)</sup>	Before first sample of the day *	Linear calibration with pH buffers 6.865 and 9.18
	Ensure sensor slope is 95% - 105% and offset < ±20 mV otherwise clean and recalibrate the sensor	
<b>Verification of measurement</b> <sup>2)</sup>	After calibration	
<b>Visual sensor verification</b>	Before Calibration. Check for blocked diaphragm and/or glass cracks	
<b>Electrolyte level maintenance</b>	Check that enough electrolyte is present before every calibration	
<b>Sensor cleaning</b>	Daily	using appropriate solution in accordance with your sample
<b>Sensor reconditioning</b>	Weekly	recondition in 0.1 M HCl
<b>Sensor storage</b>	Store in electrolyte for less than 8 hours or in storage solution for longer term storage	

System Recommendation Details			
Item	Item-No.	Pcs	Description
Seven2Go Pro S8-Biotech	30207878	1	pH/ion meter with InLab Routine Go-ISM
uGo carrying case	30122300	1	Carrying case for Seven2Go und SevenGo Duo
License EasyDirect pH	Included	1	Included Data Handling Software License
Technical Buffer pH 7.00, 30x 20 mL Sachets	51302047	2	pH buffer to cover initial 3 months demand
NIST/DIN Buffer pH 6.865, 30x 20 mL Sachets	30111137	2	pH buffer to cover initial 3 months demand
NIST/DIN Buffer pH 9.180, 30x 20 mL Sachets	30111138	2	pH buffer to cover initial 3 months demand
Storage Solution InLab, 250 mL	30111142	1	Maintenance Solution for pH sensors
StarterPac Installation	MO_StarterPac	1	Installation qualification (IQ), Initial training, Routine operation
Comprehensive Care	MO_Comp_Care	1	Preventive Maintenance on site, Remote Support, Technical Hotline, Priority Intervention, Repair Labor, Repair Parts, Travel Time

#### Remarks

1) Calibration and verification frequencies are recommended based on specified impact and measuring tolerance. Additional factors such as type of sample, number of samples and sample characteristics may result in higher or lower frequencies.

\* If the amount of samples exceeds 30 per day, it is recommended to calibrate a second time, after half the samples

2) If the verification result is outside of the specified tolerance, perform a new calibration

#### Disclaimer

The recommendations provided in this document are for information purposes only and are not binding in any way. The Recommendations are based on the information described in section "User Requirements"; other factors, which might also have an influence on the choice of systems and services, may not have been taken into account. The accuracy of measuring results is affected by the choice of meter, sensor, buffers, the properties of samples, the environment, the user, and other contributions. The recommended system and services are based on an ideal sample for the selected impact and measurement tolerance. The accuracy of your actual sample may differ from the selected tolerance and must be determined in your environment using your specific sample. The final responsibility is with the user of the system. The recommendations are provided to you without any warranties and do not extend our liability. For information on the use of your data, please visit [www.mt.com/privacypolicy](http://www.mt.com/privacypolicy).

**Date** 4/19/2023  
**Created by** parul chhaparia  
**Email** Parul.Chhaparia@mt.com